



AD Surgical

Top Surgical Essentials For Less Since 2000

January 2026

FREE SAMPLES
(for selected products)



Expertly serving *DENTAL, MEDICAL & VETERINARY PROFESSIONALS* for over 25 years

888.841.8481

www.AD-Surgical.com

CONTENTS


01	Sutures	Page 01
02	Skin Marking & Closure	Page 16
03	IV & Blood Collection	Page 17
04	Needles and Sharps	Page 19
05	Infection Control	Page 21
06	Wound Care	Page 22
07	Dental Implant Accessories	Page 23
08	Regenerative Biomaterials	Page 26
09	Dental Exam & Screening	Page 29

Retail prices listed in this catalog reflect the January 2026 printed edition only. Prices offered in our online store may vary.

NON-ABSORBABLE

Needle Size	Needle Shape	Thread Size (USP)	Thread Length	Thread Color	Item #	Qty.	Price/box
ULTRA-PREMIUM NYLON MICROSUTURES Comparable to Ethilon® * and Monosof® **							
SMALL SP-13 13mm, 3/8 Circle		7/0	18" (45cm)	Black	S-N718SP13	12 pcs	\$69.50
MICRO BV100-3 (MV-70-3) 4mm, 3/8 Circle		10/0	5" (13cm)	Black	XXS-N1005T4	12 pcs	\$225.00
MICRO BV130-4 (MV-100-4) 5mm, 3/8 Circle		9/0	5" (13cm)	Black	XXS-N905T5	12 pcs	\$225.00
MICRO BV130-5 (MV-135-5) 6.4mm, 3/8 Circle		8/0	7" (18cm)	Black	XXS-N807T6	12 pcs	\$145.00
MICRO TG-140-8 (SE-140-8) 6.4mm, 3/8 Circle		10/0	7" (18cm)	Black	XXS-N1007S6	12 pcs	\$145.00
MICRO TG-140-8 (SE-140-6) 6.4mm, 3/8 Circle		10/0	12" (30cm)	Black	XXS-N1012S6D	12 pcs	\$169.50
PREMIUM+ NYLON Comparable to Ethilon® * and Monosof® **							
LARGE PS-1 / P-14 / PC33 24mm, 3/8 Circle		2/0	27" (69cm)	Black	PLN-227R24	12 pcs	\$45.00
		3/0	27" (69cm)	Black	PLN-327R24	12 pcs	\$45.00
		4/0	27" (69cm)	Black	PLN-427R24	12 pcs	\$45.00
MEDIUM PS-2 / P-12 / PC31 19mm, 3/8 Circle		3/0	27" (69cm)	Black	PMN-327R19	12 pcs	\$45.00
		4/0	27" (69cm)	Black	PMN-427R19	12 pcs	\$45.00
		3/0	18" (45cm)	Black	PMN-318R19	12 pcs	\$39.00
		4/0	18" (45cm)	Black	PMN-418R19	12 pcs	\$39.00
MEDIUM PS-3 / P-11 / PC34 16mm, 3/8 Circle		5/0	18" (45cm)	Black	PMN-518R19	12 pcs	\$39.00
		4/0	18" (45cm)	Black	PMN-418R16	12 pcs	\$39.00
SMALL P-3 / P-13 / C3 13mm, 3/8 Circle		4/0	18" (45cm)	Black	PSN-418R13	12 pcs	\$39.00
		5/0	18" (45cm)	Black	PSN-518R13	12 pcs	\$39.00
		6/0	18" (45cm)	Black	PSN-618R13	12 pcs	\$39.00
NYLON Comparable to Ethilon® * and Monosof® **							
XX-LARGE KS / SC-2 60mm		2/0	30" (75cm)	Black	XXL-N230K60	12 pcs	\$29.00
XX-LARGE CT / GS-24 / T10 40mm, 1/2 Circle		2	30" (75cm)	Black	XXL-2N30T40	12 pcs	\$29.00



Needle Size	Needle Shape	Thread Size (USP)	Thread Length	Thread Color	Item #	Qty.	Price/box
NYLON Comparable to Ethilon® * and Monosof® **							
X-LARGE CT-1 / GS-21 / T8 36mm, 1/2 Circle		2/0	30" (75cm)	Black	XL-N230T36	12 pcs	\$29.00
X-LARGE FSL / C-16 30mm, 3/8 Circle		0	30" (75cm)	Black	XL-N030R30	12 pcs	\$29.00
LARGE FS-1 / C-14 / C7 24mm, 3/8 Circle		2/0	30" (75cm)	Black	L-N230R24	12 pcs	\$29.00
		3/0	30" (75cm)	Black	L-N330R24	12 pcs	\$29.00
		0	18" (45cm)	Black	L-N018R24	12 pcs	\$29.00
MEDIUM FS-2 / C-13 / C6 19mm, 3/8 Circle		3/0	30" (75cm)	Black	M-N330R19	12 pcs	\$29.00
		4/0	30" (75cm)	Black	M-N430R19	12 pcs	\$29.00
		3/0	18" (45cm)	Black	M-N318R19	12 pcs	\$23.00
		4/0	18" (45cm)	Black	M-N418R19	12 pcs	\$23.00
		5/0	18" (45cm)	Black	M-N518R19	12 pcs	\$23.00
		4/0	10" (25cm)	Black	M-N410R19-B	12 pcs	\$22.00
MEDIUM FS-3 / C-12 / C22 16mm, 3/8 Circle		5/0	10" (25cm)	Black	M-N510R19-B	12 pcs	\$22.00
		3/0	18" (45cm)	Black	M-N318XR19	36 pcs	\$66.00
		4/0	18" (45cm)	Black	M-N418XR19	36 pcs	\$66.00
		5/0	18" (45cm)	Black	M-N518XR19	36 pcs	\$66.00
		4/0	10" (25cm)	Black	M-N410XR19-B	36 pcs	\$63.00
SMALL P-3 / P-13 / C3 13mm, 3/8 Circle		4/0	10" (25cm)	Black	S-N410R13-B	12 pcs	\$22.00
		5/0	10" (25cm)	Black	S-N510R13-B	12 pcs	\$22.00
		6/0	10" (25cm)	Black	S-N610R13-B	12 pcs	\$22.00
		4/0	18" (45cm)	Black	S-N418R13	12 pcs	\$23.00
SMALL P-3 / P-13 / C3 13mm, 3/8 Circle		5/0	18" (45cm)	Black	S-N518R13	12 pcs	\$23.00
		6/0	18" (45cm)	Black	S-N618R13	12 pcs	\$23.00
		4/0	18" (45cm)	Black	S-N418XR13	36 pcs	\$66.00
SMALL P-3 / P-13 / C3 13mm, 3/8 Circle		5/0	18" (45cm)	Black	S-N518XR13	36 pcs	\$66.00
		6/0	18" (45cm)	Black	S-N618XR13	36 pcs	\$66.00
		6/0	18" (45cm)	Black	S-N618R11	12 pcs	\$42.50
X-SMALL P-1 / P-10 / C1 11mm, 3/8 Circle		6/0	18" (45cm)	Black	XS-N618R11	12 pcs	\$42.50

Needle Size	Needle Shape	Thread Size (USP)	Thread Length	Thread Color	Item #	Qty.	Price/box
POLYESTER Comparable to Ethibond® * and Surgidac® **							
X-LARGE FSL / C-16 30mm, 3/8 Circle		0	30" (75cm)	Green	XL-E030R30	12 pcs	\$39.50
LARGE FS / C-15 26mm, 3/8 Circle		2/0	30" (75cm)	Green	L-E230R26	12 pcs	\$39.50
LARGE FS-1 / C-14 / C7 24mm, 3/8 Circle		2/0	30" (75cm)	Green	L-E230R24	12 pcs	\$39.50
		3/0	30" (75cm)	Green	L-E330R24	12 pcs	\$39.50
MEDIUM FS-2 / C-13 / C6 19mm, 3/8 Circle		3/0	18" (45cm)	Green	M-E318R19	12 pcs	\$34.00
SMALL P-3 / P-13 / C3 13mm, 3/8 Circle		4/0	18" (45cm)	Green	S-E418R13	12 pcs	\$34.00
PREMIUM+ POLYPROPYLENE Comparable to Prolene® * and Surgipro II® **							
LARGE PS-1 / P-14 / PC33 24mm, 3/8 Circle		3/0	27" (69cm)	Blue	PLP-327R24	12 pcs	\$56.00
		4/0	27" (69cm)	Blue	PLP-427R24	12 pcs	\$56.00
MEDIUM PS-2 / P-12 / PC31 19mm, 3/8 Circle		3/0	18" (45cm)	Blue	PMP-318R19	12 pcs	\$49.00
		4/0	18" (45cm)	Blue	PMP-418R19	12 pcs	\$49.00
		5/0	18" (45cm)	Blue	PMP-518R19	12 pcs	\$49.00
MEDIUM PS-3 / P-11 / PC34 16mm, 3/8 Circle		4/0	18" (45cm)	Blue	PMP-418R16	12 pcs	\$49.00
		5/0	18" (45cm)	Blue	PMP-518R16	12 pcs	\$49.00
SMALL P-3 / P-13 / C3 13mm, 3/8 Circle		4/0	18" (45cm)	Blue	PSP-418R13	12 pcs	\$49.00
		5/0	18" (45cm)	Blue	PSP-518R13	12 pcs	\$49.00
		6/0	18" (45cm)	Blue	PSP-618R13	12 pcs	\$49.00
POLYPROPYLENE Comparable to Prolene® * and Surgipro II® **							
X-LARGE CT-1 / GS-21 / T8 36mm, 1/2 Circle		0	30" (75cm)	Blue	XL-P030T36	12 pcs	\$39.00
LARGE FS-1 / C-14 / C7 24mm, 3/8 Circle		2/0	30" (75cm)	Blue	L-P230R24	12 pcs	\$39.00

Needle Size	Needle Shape	Thread Size (USP)	Thread Length	Thread Color	Item #	Qty.	Price/box
POLYPROPYLENE Comparable to Prolene® * and Surgipro II® **							
MEDIUM FS-2 / C-13 / C6 19mm, 3/8 Circle		3/0	30" (75cm)	Blue	M-P330R19	12 pcs	\$39.00
		4/0	30" (75cm)	Blue	M-P430R19	12 pcs	\$39.00
		3/0	18" (45cm)	Blue	M-P318R19	12 pcs	\$32.00
		4/0	18" (45cm)	Blue	M-P418R19	12 pcs	\$32.00
		5/0	18" (45cm)	Blue	M-P518R19	12 pcs	\$32.00
		4/0	18" (45cm)	Blue	M-P418XR19	36 pcs	\$93.00
MEDIUM FS-3 / C-12 / C22 16mm, 3/8 Circle		4/0	18" (45cm)	Blue	M-P418R16	12 pcs	\$32.00
		5/0	18" (45cm)	Blue	M-P518R16	12 pcs	\$32.00
SMALL P-3 / P-13 / C3 13mm, 3/8 Circle		4/0	18" (45cm)	Blue	S-P418R13	12 pcs	\$32.00
		5/0	18" (45cm)	Blue	S-P518R13	12 pcs	\$32.00
		6/0	18" (45cm)	Blue	S-P618R13	12 pcs	\$32.00
		4/0	18" (45cm)	Blue	S-P418XR13	36 pcs	\$93.00
		5/0	18" (45cm)	Blue	S-P518XR13	36 pcs	\$93.00
		6/0	18" (45cm)	Blue	S-P618XR13	36 pcs	\$93.00
X-SMALL P-1 / P-10 / C1 11mm, 3/8 Circle		6/0	18" (45cm)	Blue	XS-P618R11	12 pcs	\$42.50
		7/0	18" (45cm)	Blue	XS-P718R11	12 pcs	\$42.50
PTFE Comparable to Cytoplast® *** PTFE sutures							
XX-LARGE KS / SC-2 60mm		2/0	30" (75cm)	White	XXL-F230K60	12 pcs	\$78.00
X-LARGE CT-1 / GS-21 / T8 36mm, 1/2 Circle		2/0	30" (75cm)	White	XL-F230T36	12 pcs	\$78.00
MEDIUM FS-2 / C-13 / C6 19mm, 3/8 Circle		2/0	18" (45cm)	White	M-F218R19	12 pcs	\$78.00
		3/0	18" (45cm)	White	M-F318R19	12 pcs	\$78.00
		4/0	18" (45cm)	White	M-F418R19	12 pcs	\$78.00
MEDIUM FS-3 / C-12 / C22 16mm, 3/8 Circle		3/0	18" (45cm)	White	M-F318R16	12 pcs	\$78.00
		4/0	18" (45cm)	White	M-F418R16	12 pcs	\$78.00
		5/0	18" (45cm)	White	M-F518R16	12 pcs	\$78.00
SMALL P-3 / P-13 / C3 13mm, 3/8 Circle		4/0	18" (45cm)	White	S-F418R13	12 pcs	\$78.00
		5/0	18" (45cm)	White	S-F518R13	12 pcs	\$78.00
		6/0	18" (45cm)	White	S-F618R13	12 pcs	\$78.00

Needle Size	Needle Shape	Thread Size (USP)	Thread Length	Thread Color	Item #	Qty.	Price/box
PREMIUM+ SILK Comparable to Perma-hand® * and Sofsilk® **							
MEDIUM PS-2 / P-12 / PC31 19mm, 3/8 Circle		3/0	18" (45cm)	Black	PMS-318R19	12 pcs	\$34.00
		4/0	18" (45cm)	Black	PMS-418R19	12 pcs	\$34.00
		5/0	18" (45cm)	Black	PMS-518R19	12 pcs	\$34.00
SMALL P-3 / P-13 / C3 13mm, 3/8 Circle		4/0	18" (45cm)	Black	PSS-418R13	12 pcs	\$34.00
		5/0	18" (45cm)	Black	PSS-518R13	12 pcs	\$34.00
		6/0	18" (45cm)	Black	PSS-618R13	12 pcs	\$34.00
SILK Comparable to Perma-hand® * and Sofsilk® **							
XX-LARGE KS / SC-2 60mm		2/0	30" (75cm)	Black	XXL-S230K60	12 pcs	\$26.00
MEDIUM FS-2 / C-13 / C6 19mm, 3/8 Circle		3/0	30" (75cm)	Black	M-S330R19	12 pcs	\$26.00
		4/0	30" (75cm)	Black	M-S430R19	12 pcs	\$26.00
		3/0	18" (45cm)	Black	M-S318R19	12 pcs	\$22.00
MEDIUM FS-3 / C-12 / C22 16mm, 3/8 Circle		4/0	18" (45cm)	Black	M-S418R19	12 pcs	\$22.00
		5/0	18" (45cm)	Black	M-S518R19	12 pcs	\$22.00
		4/0	18" (45cm)	Black	M-S418R16	12 pcs	\$22.00
SMALL P-3 / P-13 / C3 13mm, 3/8 Circle		5/0	18" (45cm)	Black	M-S518R16	12 pcs	\$22.00
		4/0	18" (45cm)	Black	S-S418R13	12 pcs	\$22.00
		5/0	18" (45cm)	Black	S-S518R13	12 pcs	\$22.00
6/0	18" (45cm)	Black	S-S618R13	12 pcs	\$22.00		
PREMIUM+ PDO Comparable to PDS II® * and Maxon® **							
LARGE PS-1 / P-14 / PC33 24mm, 3/8 Circle		2/0	27" (69cm)	Violet	PLD-227R24	12 pcs	\$59.50
		3/0	27" (69cm)	Violet	PLD-327R24	12 pcs	\$59.50
		4/0	27" (69cm)	Violet	PLD-427R24	12 pcs	\$59.50
MEDIUM PS-2 / P-12 / PC31 19mm, 3/8 Circle		2/0	27" (69cm)	Violet	PMD-227R19	12 pcs	\$59.50
		3/0	27" (69cm)	Violet	PMD-327R19	12 pcs	\$59.50
		4/0	27" (69cm)	Violet	PMD-427R19	12 pcs	\$59.50
		2/0	18" (45cm)	Violet	PMD-218R19	12 pcs	\$52.00
		3/0	18" (45cm)	Violet	PMD-318R19	12 pcs	\$52.00
		4/0	18" (45cm)	Violet	PMD-418R19	12 pcs	\$52.00
MEDIUM PS-3 / P-11 / PC34 16mm, 3/8 Circle		5/0	18" (45cm)	Violet	PMD-518R19	12 pcs	\$52.00
		4/0	18" (45cm)	Violet	PMD-418R16	12 pcs	\$52.00
		4/0	18" (45cm)	Violet	PSD-418R13	12 pcs	\$52.00
SMALL P-3 / P-13 / C3		5/0	18" (45cm)	Violet	PSD-518R13	12 pcs	\$52.00

Needle Size	Needle Shape	Thread Size (USP)	Thread Length	Thread Color	Item #	Qty.	Price/box
PDO Comparable to PDS II® * and Maxon® **							
XX-LARGE KS / SC-2 60mm Straight		4/0	30" (75cm)	Violet	XXL-D430K60	12 pcs	\$52.00
XX-LARGE CP / GS-12 / C14 40mm, 1/2 Circle		2	30" (75cm)	Violet	XXL-2D30R40	12 pcs	\$52.00
		1	30" (75cm)	Violet	XXL-1D30R40	12 pcs	\$52.00
X-LARGE CT-1 / GS-21 / T8 36mm, 1/2 Circle		0	30" (75cm)	Violet	XL-D030T36	12 pcs	\$52.00
		2/0	30" (75cm)	Violet	XL-D230T36	12 pcs	\$52.00
X-LARGE FSL / C-16 30mm, 3/8 Circle		0	30" (75cm)	Violet	XL-D030R30	12 pcs	\$52.00
LARGE SH (V-20) 26mm, 1/2 Circle		2/0	30" (75cm)	Violet	L-D230T26	12 pcs	\$48.00
		3/0	30" (75cm)	Violet	L-D330T26	12 pcs	\$48.00
LARGE FS-1 / C-14 / C7 24mm, 3/8 Circle		2/0	30" (75cm)	Violet	L-D230R24	12 pcs	\$48.00
		3/0	30" (75cm)	Violet	L-D330R24	12 pcs	\$48.00
		2/0	30" (75cm)	Violet	L-D230XR24	36 pcs	\$136.00
		3/0	30" (75cm)	Violet	L-D330XR24	36 pcs	\$136.00
MEDIUM FS-2 / C-13 / C6 19mm, 3/8 Circle		2/0	30" (75cm)	Violet	M-D230R19	12 pcs	\$48.00
		3/0	30" (75cm)	Violet	M-D330R19	12 pcs	\$48.00
		4/0	30" (75cm)	Violet	M-D430R19	12 pcs	\$48.00
		3/0	18" (45cm)	Violet	M-D318R19	12 pcs	\$44.00
		4/0	18" (45cm)	Violet	M-D418R19	12 pcs	\$44.00
MEDIUM RB-1 / CV-23 / T28 17mm, 1/2 Circle		5/0	18" (45cm)	Violet	M-D518R19	12 pcs	\$44.00
		3/0	30" (75cm)	Violet	M-D330T17	12 pcs	\$48.00
		3/0	30" (75cm)	Violet	M-D330T17	12 pcs	\$48.00
MEDIUM FS-3 / C-12 / C22 16mm, 3/8 Circle		4/0	18" (45cm)	Violet	M-D418R16	12 pcs	\$44.00
		5/0	18" (45cm)	Violet	M-D518R16	12 pcs	\$44.00
SMALL P-3 / P-13 / C3 13mm, 3/8 Circle		4/0	18" (45cm)	Violet	S-D418R13	12 pcs	\$44.00
		5/0	18" (45cm)	Violet	S-D518R13	12 pcs	\$44.00
		6/0	18" (45cm)	Violet	S-D618R13	12 pcs	\$44.00
Knotless Barbed PDO Comparable to Stratafix PDO® * and V-Loc® **							
LARGE SH (V-20) 26mm, 1/2 Circle		2/0	18" (45cm)	Violet	B-D218T26	6 pcs	\$139.50
MEDIUM PS-2 / P-12 / PC31 19mm, 3/8 Circle		3/0	12" (30cm)	Violet	B-PD312R19	6 pcs	\$139.50

Needle Size	Needle Shape	Thread Size (USP)	Thread Length	Thread Color	Item #	Qty.	Price/box
ULTRA-PREMIUM PGA MICROSUTURES Comparable to Vicryl® * and Polysorb® **							
MICRO TG-140-8 (SE-140-6) 6.4mm, 3/8 Circle		6/0	12" (30cm)	Violet	XXS-G612S6D	12 pcs	\$169.50
		8/0	12" (30cm)	Violet	XXS-G812S6D	12 pcs	\$169.50
PREMIUM+ PGA Comparable to Vicryl® * and Polysorb® **							
XX-LARGE CT / GS-24 / T10 40mm, 1/2 Circle		2/0	27" (69cm)	Violet	PXXL-G227T40	12 pcs	\$59.50
X-LARGE CT-1 / GS-21 / T8 36mm, 1/2 Circle		2/0	27" (69cm)	Undyed	PXLGU-227T36	12 pcs	\$59.50
LARGE SH (V-20) 26mm, 1/2 Circle		3/0	27" (69cm)	Undyed	PLGU-327T26	12 pcs	\$59.50
LARGE PS-1 / P-14 / PC33 24mm, 3/8 Circle		2/0	27" (69cm)	Undyed	PLGU-227R24	12 pcs	\$59.50
		3/0	27" (69cm)	Undyed	PLGU-327R24	12 pcs	\$59.50
		3/0	27" (69cm)	Violet	PLG-327R24	12 pcs	\$59.50
		4/0	27" (69cm)	Violet	PLG-427R24	12 pcs	\$59.50
MEDIUM PS-2 / P-12 / PC31 19mm, 3/8 Circle		3/0	27" (69cm)	Violet	PMG-327R19	12 pcs	\$59.50
		4/0	27" (69cm)	Violet	PMG-427R19	12 pcs	\$59.50
		3/0	27" (69cm)	Undyed	PMGU-327R19	12 pcs	\$59.50
		4/0	27" (69cm)	Undyed	PMGU-427R19	12 pcs	\$59.50
		3/0	18" (45cm)	Violet	PMG-318R19	12 pcs	\$52.00
		4/0	18" (45cm)	Violet	PMG-418R19	12 pcs	\$52.00
		5/0	18" (45cm)	Violet	PMG-518R19	12 pcs	\$52.00
		3/0	18" (45cm)	Undyed	PMGU-318R19	12 pcs	\$52.00
MEDIUM PS-3 / P-11 / PC34 16mm, 3/8 Circle		4/0	18" (45cm)	Undyed	PMGU-418R16	12 pcs	\$52.00
		5/0	18" (45cm)	Undyed	PMGU-518R16	12 pcs	\$52.00
		4/0	18" (45cm)	Violet	PMG-418R16	12 pcs	\$52.00
SMALL P-3 / P-13 / C3 13mm, 3/8 Circle		5/0	18" (45cm)	Violet	PMG-518R16	12 pcs	\$52.00
		4/0	18" (45cm)	Undyed	PSGU-418R13	12 pcs	\$52.00
		5/0	18" (45cm)	Undyed	PSGU-518R13	12 pcs	\$52.00
		6/0	18" (45cm)	Undyed	PSGU-618R13	12 pcs	\$52.00
		4/0	18" (45cm)	Violet	PSG-418R13	12 pcs	\$52.00
		5/0	18" (45cm)	Violet	PSG-518R13	12 pcs	\$52.00
6/0	18" (45cm)	Violet	PSG-618R13	12 pcs	\$52.00		

Needle Size	Needle Shape	Thread Size (USP)	Thread Length	Thread Color	Item #	Qty.	Price/box
PGA Comparable to Vicryl® *and Polysorb® **							
XX-LARGE CT / GS-24 / T10 40mm, 1/2 Circle		1	30" (75cm)	Violet	XXL-1G30T40	12 pcs	\$52.00
		0	30" (75cm)	Violet	XXL-G030T40	12 pcs	\$52.00
		2/0	30" (75cm)	Violet	XXL-G230T40	12 pcs	\$52.00
		1	30" (75cm)	Undyed	XXL-1G30T40-U	12 pcs	\$52.00
		0	30" (75cm)	Undyed	XXL-G030T40-U	12 pcs	\$52.00
		2/0	30" (75cm)	Undyed	XXL-G230T40-U	12 pcs	\$52.00
X-LARGE CT-1 / GS-21 / T8 36mm, 1/2 Circle		0	30" (75cm)	Violet	XL-G030T36	12 pcs	\$52.00
		2/0	30" (75cm)	Undyed	XL-G230T36-U	12 pcs	\$52.00
X-LARGE FSL / C-16 30mm, 3/8 Circle		0	30" (75cm)	Violet	XL-G030R30	12 pcs	\$52.00
LARGE SH (V-20) 26mm, 1/2 Circle		2/0	30" (75cm)	Violet	L-G230T26	12 pcs	\$48.00
		2/0	30" (75cm)	Undyed	L-G230T26-U	12 pcs	\$48.00
		3/0	30" (75cm)	Violet	L-G330T26	12 pcs	\$48.00
LARGE FS-1 / C-14 / C7 24mm, 3/8 Circle		2/0	30" (75cm)	Violet	L-G230R24	12 pcs	\$48.00
		3/0	30" (75cm)	Violet	L-G330R24	12 pcs	\$48.00
		2/0	30" (75cm)	Undyed	L-G230R24-U	12 pcs	\$48.00
		3/0	30" (75cm)	Undyed	L-G330R24-U	12 pcs	\$48.00
				3/0	30" (75cm)	Violet	M-G330R19
		4/0	30" (75cm)	Violet	M-G430R19	12 pcs	\$48.00
		3/0	30" (75cm)	Undyed	M-G330R19-U	12 pcs	\$48.00
		4/0	30" (75cm)	Undyed	M-G430R19-U	12 pcs	\$48.00
MEDIUM FS-2 / C-13 / C6 19mm, 3/8 Circle		3/0	18" (45cm)	Violet	M-G318R19	12 pcs	\$44.00
		4/0	18" (45cm)	Violet	M-G418R19	12 pcs	\$44.00
		5/0	18" (45cm)	Violet	M-G518R19	12 pcs	\$44.00
		3/0	18" (45cm)	Undyed	M-G318R19-U	12 pcs	\$44.00
		4/0	18" (45cm)	Undyed	M-G418R19-U	12 pcs	\$44.00
		5/0	18" (45cm)	Undyed	M-G518R19-U	12 pcs	\$44.00
		3/0	18" (45cm)	Violet	M-G318XR19	36 pcs	\$125.00
		4/0	18" (45cm)	Violet	M-G418XR19	36 pcs	\$125.00
		3/0	18" (45cm)	Undyed	M-G318XR19-U	36 pcs	\$125.00
		4/0	18" (45cm)	Undyed	M-G418XR19-U	36 pcs	\$125.00
MEDIUM RB-1 / CV-23 / T28 17mm, 1/2 Circle		3/0	30" (75cm)	Undyed	M-G330T17-U	12 pcs	\$48.00
		4/0	30" (75cm)	Violet	M-G430T17	12 pcs	\$48.00
		4/0	30" (75cm)	Undyed	M-G430T17-U	12 pcs	\$48.00

Needle Size	Needle Shape	Thread Size (USP)	Thread Length	Thread Color	Item #	Qty.	Price/box
PGA Comparable to Vicryl® *and Polysorb® **							
MEDIUM FS-3 / C-12 / C22 16mm, 3/8 Circle		4/0	18" (45cm)	Violet	M-G418R16	12 pcs	\$44.00
		5/0	18" (45cm)	Violet	M-G518R16	12 pcs	\$44.00
		4/0	18" (45cm)	Undyed	M-G418R16-U	12 pcs	\$44.00
		5/0	18" (45cm)	Undyed	M-G518R16-U	12 pcs	\$44.00
MEDIUM PC-3 / PC-11 16mm, 3/8 Circle		4/0	18" (45cm)	Undyed	M-G418A16-U	12 pcs	\$44.00
		4/0	18" (45cm)	Violet	S-G418R13	12 pcs	\$44.00
		5/0	18" (45cm)	Violet	S-G518R13	12 pcs	\$44.00
		6/0	18" (45cm)	Violet	S-G618R13	12 pcs	\$44.00
SMALL P-3 / P-13 / C3 13mm, 3/8 Circle		4/0	18" (45cm)	Undyed	S-G418R13-U	12 pcs	\$44.00
		5/0	18" (45cm)	Undyed	S-G518R13-U	12 pcs	\$44.00
		6/0	18" (45cm)	Undyed	S-G618R13-U	12 pcs	\$44.00
		4/0	18" (45cm)	Violet	S-G418XR13	36 pcs	\$125.00
		4/0	18" (45cm)	Undyed	S-G418XR13-U	36 pcs	\$125.00
		5/0	18" (45cm)	Undyed	S-G518XR13-U	36 pcs	\$125.00
X-SMALL P-1 / P-10 / C1 11mm, 3/8 Circle		5/0	18" (45cm)	Undyed	XS-G518R11-U	12 pcs	\$54.00
		6/0	18" (45cm)	Undyed	XS-G618R11-U	12 pcs	\$54.00
		7/0	18" (45cm)	Violet	XS-G718R11	12 pcs	\$54.00
PREMIUM+ RAPID PGA Comparable to Vicryl Rapide* and Velosorb Fast**							
MEDIUM PS-2 / P-12 / PC31 19mm, 3/8 Circle		4/0	18" (45cm)	Undyed	PMK-418R19	12 pcs	\$54.00
SMALL P-3 / P-13 / C3 13mm, 3/8 Circle		4/0	18" (45cm)	Undyed	PSK-418R13	12 pcs	\$54.00
		5/0	18" (45cm)	Undyed	PSK-518R13	12 pcs	\$54.00
RAPID PGA Comparable to Vicryl Rapide* and Velosorb Fast**							
MEDIUM FS-2 / C-13 / C6 19mm, 3/8 Circle		4/0	30" (75cm)	Undyed	M-K430R19	12 pcs	\$52.00
		3/0	18" (45cm)	Undyed	M-K318R19	12 pcs	\$44.00
		4/0	18" (45cm)	Undyed	M-K418R19	12 pcs	\$44.00
		5/0	18" (45cm)	Undyed	M-K518R19	12 pcs	\$44.00
SMALL P-3 / P-13 / C3 13mm, 3/8 Circle		4/0	18" (45cm)	Undyed	S-K418R13	12 pcs	\$44.00
		5/0	18" (45cm)	Undyed	S-K518R13	12 pcs	\$44.00
X-SMALL P-1 / P-10 / C1 11mm, 3/8 Circle		5/0	18" (45cm)	Undyed	XS-K518R11	12 pcs	\$54.00

Needle Size	Needle Shape	Thread Size (USP)	Thread Length	Thread Color	Item #	Qty.	Price/box
PREMIUM+ PGCL Comparable to Monocryl® * and Biosyn®**							
LARGE PS-1 / P-14 / PC33 24mm, 3/8 Circle		3/0	27" (69cm)	Violet	PLQC-327R24	12 pcs	\$59.50
		4/0	27" (69cm)	Violet	PLQC-427R24	12 pcs	\$59.50
		2/0	27" (69cm)	Undyed	PLQ-227R24	12 pcs	\$59.50
		3/0	27" (69cm)	Undyed	PLQ-327R24	12 pcs	\$59.50
		4/0	27" (69cm)	Undyed	PLQ-427R24	12 pcs	\$59.50
MEDIUM PS-2 / P-12 / PC31 19mm, 3/8 Circle		3/0	27" (69cm)	Undyed	PMQ-327R19	12 pcs	\$59.50
		4/0	27" (69cm)	Undyed	PMQ-427R19	12 pcs	\$59.50
		3/0	18" (45cm)	Violet	PMQC-318R19	12 pcs	\$52.00
		4/0	18" (45cm)	Violet	PMQC-418R19	12 pcs	\$52.00
		5/0	18" (45cm)	Violet	PMQC-518R19	12 pcs	\$52.00
SMALL P-3 / P-13 / C3 13mm, 3/8 Circle		3/0	18" (45cm)	Undyed	PMQ-318R19	12 pcs	\$52.00
		4/0	18" (45cm)	Undyed	PMQ-418R19	12 pcs	\$52.00
		5/0	18" (45cm)	Undyed	PMQ-518R19	12 pcs	\$52.00
		4/0	18" (45cm)	Violet	PSQC-418R13	12 pcs	\$52.00
		5/0	18" (45cm)	Violet	PSQC-518R13	12 pcs	\$52.00
SMALL TF / CV-22 / T22 13mm, 1/2 Circle		4/0	18" (45cm)	Undyed	PSQ-418R13	12 pcs	\$52.00
		5/0	18" (45cm)	Undyed	PSQ-518R13	12 pcs	\$52.00
		6/0	18" (45cm)	Undyed	PSQ-618R13	12 pcs	\$52.00
5/0	18" (45cm)	Undyed	PSQ-518T13	12 pcs	\$52.00		

Needle Size	Needle Shape	Thread Size (USP)	Thread Length	Thread Color	Item #	Qty.	Price/box
PGCL Comparable to Monocryl® * and Biosyn®**							
X-LARGE CT-1 / GS-21 / T8 36mm, 1/2 Circle		2/0	30" (75cm)	Undyed	XL-Q230T36	12 pcs	\$52.00
		3/0	30" (75cm)	Undyed	M-Q330R19	12 pcs	\$52.00
LARGE FS-1 / C-14 / C7 24mm, 3/8 Circle		3/0	30" (75cm)	Undyed	L-Q230R24	12 pcs	\$52.00
		3/0	30" (75cm)	Undyed	L-Q330R24	12 pcs	\$52.00
		4/0	30" (75cm)	Undyed	L-Q430R24	12 pcs	\$52.00
		4/0	30" (75cm)	Undyed	M-Q430R19	12 pcs	\$52.00
		3/0	18" (45cm)	Undyed	M-Q318R19	12 pcs	\$44.00
MEDIUM FS-2 / C-13 / C6 19mm, 3/8 Circle		4/0	18" (45cm)	Undyed	M-Q418R19	12 pcs	\$44.00
		5/0	18" (45cm)	Undyed	M-Q518R19	12 pcs	\$44.00
		3/0	18" (45cm)	Violet	M-Q318R19-C	12 pcs	\$44.00
		4/0	18" (45cm)	Violet	M-Q418R19-C	12 pcs	\$44.00
		5/0	18" (45cm)	Violet	M-Q518R19-C	12 pcs	\$44.00
SMALL P-3 / P-13 / C3 13mm, 3/8 Circle		3/0	30" (75cm)	Undyed	M-Q330XR19	36 pcs	\$150.00
		4/0	30" (75cm)	Undyed	M-Q430XR19	36 pcs	\$150.00
		3/0	18" (45cm)	Undyed	M-Q318XR19	36 pcs	\$125.00
		4/0	18" (45cm)	Undyed	M-Q418XR19	36 pcs	\$125.00
		4/0	18" (45cm)	Violet	M-Q418XR19-C	36 pcs	\$125.00
		4/0	18" (45cm)	Undyed	S-Q418R13	12 pcs	\$44.00
		5/0	18" (45cm)	Undyed	S-Q518R13	12 pcs	\$44.00
		6/0	18" (45cm)	Undyed	S-Q618R13	12 pcs	\$44.00
		4/0	18" (45cm)	Violet	S-Q418R13-C	12 pcs	\$44.00
		5/0	18" (45cm)	Violet	S-Q518R13-C	12 pcs	\$44.00
6/0	18" (45cm)	Violet	S-Q618R13-C	12 pcs	\$44.00		
X-SMALL P-1 / P-10 / C1 11mm, 3/8 Circle		4/0	18" (45cm)	Undyed	S-Q418XR13	36 pcs	\$125.00
		5/0	18" (45cm)	Undyed	S-Q518XR13	36 pcs	\$125.00
		5/0	18" (45cm)	Violet	S-Q518XR13-C	36 pcs	\$125.00
5/0	18" (45cm)	Violet	XS-Q518R11-C	12 pcs	\$54.00		

* Registered trademarks of Ethicon

** Registered trademarks of Covidien

*** Registered trademark of Osteogenics Biomedical

Needle Size	Needle Shape	Thread Size (USP)	Thread Length	Thread Color	Item #	Qty.	Price/box
PREMIUM+ FAST GUT							
SMALL P-3 / P-13 / C3 13mm, 3/8 Circle		5/0	18" (45cm)	Beige	PSRJ-518R13	12 pcs	\$56.00
		6/0	18" (45cm)	Beige	PSRJ-618R13	12 pcs	\$56.00
SMALL PC-1 13mm, 3/8 Circle		5/0	18" (45cm)	Beige	PSRJ-518A13	12 pcs	\$56.00
PREMIUM+ PLAIN GUT							
LARGE PS-1 / P-14 / PC33 24mm, 3/8 Circle		3/0	27" (69cm)	Beige	PLJ-327R24	12 pcs	\$49.50
		4/0	27" (69cm)	Beige	PMJ-427R19	12 pcs	\$49.50
MEDIUM PS-2 / P-12 / PC31 19mm, 3/8 Circle		3/0	27" (69cm)	Beige	PMJ-327R19	12 pcs	\$49.50
		4/0	27" (69cm)	Beige	PMJ-427R19	12 pcs	\$49.50
		3/0	18" (45cm)	Beige	PMJ-318R19	12 pcs	\$48.00
MEDIUM PS-3 / P-11 / PC34 16mm, 3/8 Circle		4/0	18" (45cm)	Beige	PMJ-418R19	12 pcs	\$48.00
		5/0	18" (45cm)	Beige	PMJ-518R16	12 pcs	\$48.00
SMALL P-3 / P-13 / C3 13mm, 3/8 Circle		4/0	18" (45cm)	Beige	PSJ-418R13	12 pcs	\$48.00
		5/0	18" (45cm)	Beige	PSJ-518R13	12 pcs	\$48.00
		6/0	18" (45cm)	Beige	PSJ-618R13	12 pcs	\$48.00
PLAIN GUT							
LARGE FS-1 / C-14 / C7 24mm, 3/8 Circle		2/0	30" (75cm)	Beige	L-J230R24	12 pcs	\$42.00
		3/0	30" (75cm)	Beige	M-J330R19	12 pcs	\$42.00
MEDIUM FS-2 / C-13 / C6 19mm, 3/8 Circle		4/0	30" (75cm)	Beige	M-J430R19	12 pcs	\$42.00
		3/0	18" (45cm)	Beige	M-J318R19	12 pcs	\$39.50
		4/0	18" (45cm)	Beige	M-J418R19	12 pcs	\$39.50
		5/0	18" (45cm)	Beige	M-J518R19	12 pcs	\$39.50
MEDIUM FS-3 / C-12 / C22 16mm, 3/8 Circle		4/0	18" (45cm)	Beige	M-J418R16	12 pcs	\$39.50
		5/0	18" (45cm)	Beige	M-J518R16	12 pcs	\$39.50
SMALL P-3 / P-13 / C3 13mm, 3/8 Circle		4/0	18" (45cm)	Beige	S-J418R13	12 pcs	\$39.50
		5/0	18" (45cm)	Beige	S-J518R13	12 pcs	\$39.50
		6/0	18" (45cm)	Beige	S-J618R13	12 pcs	\$39.50

Needle Size	Needle Shape	Thread Size (USP)	Thread Length	Thread Color	Item #	Qty.	Price/box
PREMIUM+ CHROMIC GUT							
LARGE PS-1 / P-14 / PC33 24mm, 3/8 Circle		2/0	27" (69cm)	Brown	PLC-227R24	12 pcs	\$49.50
		3/0	27" (69cm)	Brown	PLC-327R24	12 pcs	\$49.50
MEDIUM PS-2 / P-12 / PC31 19mm, 3/8 Circle		3/0	27" (69cm)	Brown	PMC-327R19	12 pcs	\$49.50
		4/0	27" (69cm)	Brown	PMC-427R19	12 pcs	\$49.50
		3/0	18" (45cm)	Brown	PMC-318R19	12 pcs	\$48.00
MEDIUM PS-3 / P-11 / PC34 16mm, 3/8 Circle		4/0	18" (45cm)	Brown	PMC-418R19	12 pcs	\$48.00
		5/0	18" (45cm)	Brown	PMC-518R16	12 pcs	\$48.00
SMALL P-3 / P-13 / C3 13mm, 3/8 Circle		4/0	18" (45cm)	Brown	PSC-418R13	12 pcs	\$48.00
		5/0	18" (45cm)	Brown	PSC-518R13	12 pcs	\$48.00
		6/0	18" (45cm)	Brown	PSC-618R13	12 pcs	\$48.00
CHROMIC GUT							
LARGE FS-1 / C-14 / C7 24mm, 3/8 Circle		2/0	30" (75cm)	Brown	L-C230R24	12 pcs	\$42.00
		3/0	30" (75cm)	Brown	M-C330R19	12 pcs	\$42.00
MEDIUM FS-2 / C-13 / C6 19mm, 3/8 Circle		4/0	30" (75cm)	Brown	M-C430R19	12 pcs	\$42.00
		3/0	18" (45cm)	Brown	M-C318R19	12 pcs	\$39.50
		4/0	18" (45cm)	Brown	M-C418R19	12 pcs	\$39.50
		5/0	18" (45cm)	Brown	M-C518R19	12 pcs	\$39.50
		4/0	18" (45cm)	Brown	M-C418R16	12 pcs	\$39.50
MEDIUM FS-3 / C-12 / C22 16mm, 3/8 Circle		5/0	18" (45cm)	Brown	M-C518R16	12 pcs	\$39.50
		4/0	18" (45cm)	Brown	S-C418R13	12 pcs	\$39.50
SMALL P-3 / P-13 / C3 13mm, 3/8 Circle		5/0	18" (45cm)	Brown	S-C518R13	12 pcs	\$39.50
		6/0	18" (45cm)	Brown	S-C618R13	12 pcs	\$39.50

UNIFY® VET SUTURES For Veterinary Use Only

Needle Size	Needle Shape	Thread Size (USP)	Thread Length	Thread Color	Item #	Qty.	Price/box
VET NYLON Comparable to Ethilon® * and Monosof® **							
X-LARGE CP-1 / GS-11 36mm, 1/2 Circle		0	30" (75cm)	Black	VETN-030R36	12 pcs	\$20.00
X-LARGE FSL / C-16 30mm, 3/8 Circle		0	30" (75cm)	Black	VETN-030R30	12 pcs	\$20.00
LARGE FS-1 / C-14 / C7 24mm, 3/8 Circle		2/0 3/0	30" (75cm) 30" (75cm)	Black Black	VETN-230R24 VETN-330R24	12 pcs 12 pcs	\$20.00 \$20.00
MEDIUM FS-2 / C-13 / C6 19mm, 3/8 Circle		4/0	30" (75cm)	Black	VETN-430R19	12 pcs	\$20.00
VET PGA Comparable to Vicryl® * and Polysorb® **							
X-LARGE CP-1 / GS-11 36mm, 1/2 Circle		0	30" (75cm)	Violet	VETG-030R36	12 pcs	\$34.00
X-LARGE FSL / C-16 30mm, 3/8 Circle		0	30" (75cm)	Violet	VETG-030R30	12 pcs	\$34.00
LARGE FS-1 / C-14 / C7 24mm, 3/8 Circle		2/0 3/0	30" (75cm) 30" (75cm)	Violet Violet	VETG-230R24 VETG-330R24	12 pcs 12 pcs	\$34.00 \$34.00
MEDIUM FS-2 / C-13 / C6 19mm, 3/8 Circle		4/0	30" (75cm)	Violet	VETG-430R19	12 pcs	\$34.00
VET PDO Comparable to PDS II® * and Maxon® **							
X-LARGE CP-1 / GS-11 36mm, 1/2 Circle		0	30" (75cm)	Violet	VETD-030R36	12 pcs	\$38.00
X-LARGE FSL / C-16 30mm, 3/8 Circle		0	30" (75cm)	Violet	VETD-030R30	12 pcs	\$38.00

UNIFY® VET SUTURES For Veterinary Use Only

Needle Size	Needle Shape	Thread Size (USP)	Thread Length	Thread Color	Item #	Qty.	Price/box
VET PDO Comparable to PDS II® * and Maxon® **							
LARGE FS-1 / C-14 / C7 24mm, 3/8 Circle		2/0 3/0	30" (75cm) 30" (75cm)	Violet Violet	VETD-230R24 VETD-330R24	12 pcs 12 pcs	\$38.00 \$38.00
MEDIUM FS-2 / C-13 / C6 19mm, 3/8 Circle		4/0	30" (75cm)	Violet	VETD-430R19	12 pcs	\$38.00
VET PGCL Comparable to Monocryl® * and Biosyn®**							
X-LARGE CP-1 / GS-11 36mm, 1/2 Circle		0	30" (75cm)	Violet	VETQ-030R36-C	12 pcs	\$38.00
X-LARGE FSL / C-16 30mm, 3/8 Circle		0	30" (75cm)	Violet	VETQ-030R30-C	12 pcs	\$38.00
LARGE FS-1 / C-14 / C7 24mm, 3/8 Circle		2/0 3/0	30" (75cm) 30" (75cm)	Violet Violet	VETQ-230R24-C VETQ-330R24-C	12 pcs 12 pcs	\$38.00 \$38.00
MEDIUM FS-2 / C-13 / C6 19mm, 3/8 Circle		4/0	30" (75cm)	Violet	VETQ-430R19-C	12 pcs	\$38.00
VET Chromic Gut							
X-LARGE CP-1 / GS-11 36mm, 1/2 Circle		0	30" (75cm)	Brown	VETC-030R36	12 pcs	\$34.00
X-LARGE FSL / C-16 30mm, 3/8 Circle		0	30" (75cm)	Brown	VETC-030R30	12 pcs	\$34.00
LARGE FS-1 / C-14 / C7 24mm, 3/8 Circle		2/0 3/0	30" (75cm) 30" (75cm)	Brown Brown	VETC-230R24 VETC-330R24	12 pcs 12 pcs	\$34.00 \$34.00
MEDIUM FS-2 / C-13 / C6 19mm, 3/8 Circle		4/0	30" (75cm)	Brown	VETC-430R19	12 pcs	\$34.00

Surgical Skin Markers (Violet Ink)



AD SURGICAL Surgical Skin Marker is designed to mark surgical incisions and anatomical sites on a patient's skin prior to surgical procedure to promote correct site surgery. Each sterile and ready-to-use surgical skin marker includes a 6" ruler and 6 blank labels. Two tip sizes available for different user and/or surgical procedure needs. Natural gentian violet plant extract ink is hypoallergenic and non-toxic. For single use only. Sterilized by gamma radiation.

Size	Item#	BOX OF 50
0.5mm Fine Tip	SKM-0005	\$49.50
1.0mm Regular Tip	SKM-1000	\$49.50
0.5mm + 1.0mm Tips	SKM-0510	\$49.50

Aesthetic Skin Markers (Erasable White Ink)



AD SURGICAL Aesthetic Skin Marker has erasable white ink that can be easily wiped off without washing or scrubbing. Ideal for consultation appointments and/or non-surgical aesthetic procedures like Botox®, dermal fillers, Kybella® injection sites, and ultrasonic liposuction. Available only in 1.0mm (Regular) tip. For single use only. Sterilized by gamma radiation.

Size	Item#	BOX OF 50
1.0mm Tip	SKM-WH1000	\$49.50

Topical Skin Adhesives



GLUSEAL® Topical Skin Adhesives are liquid adhesive bandages intended to cover cuts, scrapes, burns and minor irritations of the skin and protect them from infection. Each 5ml multi-use kit contains enough tissue adhesive for up to 50 applications. Available in multi-use and single-use versions.

Item # GLU-550
Packing: One Multi-Use Kit (1 x 5ml skin adhesive, 50 x micro-pipettes, 1 x autoclavable tray, 25 x foam tipped applicators)

Price: \$130.00 /kit



Item # GLU-55012
Packing: Single-use kit
 12 x 0.2ml /box
Price: \$130 /box

Disposable Skin Staplers



EASYTAP® Disposable Skin Staplers are lever action staplers with an ergonomic design and strong construction to provide reliable performance during surgery. Each stapler is preloaded with 35 wide (6.5mm x 4.7mm) staples. Individually sealed & sterilized by EO Gas.

Item # SM-35W
Packing: 6 sealed staplers per box
Price: \$79.50 /box

Suturing Practice Skin Pads

AD SURGICAL Suturing Practice Skin Pads are made from a durable high-grade silicone material that can be sutured repeatedly for practice. Each ultra-realistic pad mimics closely the 3 main layers of the skin (skin/ fat/ muscle) in texture and appearance. Available with or without pre-made wounds. Made in USA.

5" x 7" pad with wounds (no case)



Item # T57-3SWP
Packing: 1 PC
Price: \$49.00 each

5" x 5" pad without wound (with case)



Item # T57-Skin
Packing: 1 PC
Price: \$59.00 each

Cassette Sutures (VET use only)



UNIFY® PDO Cassette Sutures are high-quality suture thread materials provided in a convenient and economical delivery system. All cassette sutures are intended for veterinary use only. Sterilized by EO Gas.

PDO				
0	15 meters	Violet	VCD-L0-15	\$54.00
2/0	15 meters	Violet	VCD-20-15	\$54.00
3/0	15 meters	Violet	VCD-30-15	\$54.00

Blood Collection Sets (with Holder)



KBMED Safety Blood Collection Set with Holder features wide, flexible wings that provide greater control while positioning the needle. The safety design helps prevent needlestick injuries using a safety shield that pushes over the needle following venipuncture. Pre-attached holder meets OSHA single-use holder standard. Sterilized by EO gas.

Size:	Item#
21G x 3/4"	#SBCN21-1H
23G x 3/4"	#SBCN23K-1H
25G x 3/4"	#SBCN25-1H

Packing: Box of 25
Price: \$35.00 /box

Blood Collection Sets (w/o Holder)



KB MED Blood Collection Sets (without holder) allow users to collect blood draw samples with minimal risk due to their safety mechanism. They feature wide, flexible wings that provide greater control while positioning the needle, making one-handed application easier. The safety design helps prevent needlestick injuries using a safety shield that pushes over the needle following venipuncture. Sterilized by EO gas.

Size:	Item#
21G x 3/4"	#SBCN21-1
23G x 3/4"	#SBCN23-1
25G x 3/4"	#SBCN25-1

Packing: Box of 50
Price: \$45.00 /box

Winged IV Catheters (w/o safety)



SHARD Premium+ Winged IV Catheters (without safety) are designed to perform in delicate vasculature, delivering maximum efficiency and patient comfort. The smooth catheter-to-needle transition results in atraumatic insertion for the patient. Catheter features a flashback window for quick visualization of blood. Hub is color-coded for rapid identification of catheter size, and the winged design offers enhanced control and stability. Sterilized by EO Gas.

Size:	Item#
18G x 1.25" (32mm)	#CAT-1832
20G x 1" (25mm)	#CAT-2025
20G x 1.25" (32mm)	#CAT-2032
22G x 1" (25mm)	#CAT-2225
24G x 0.75" (19mm)	#CAT-2419

Packing: Box of 50
Price: \$59.50 /box

Prefilled Saline IV Flush Syringes



SAFE-T-FLUSH Prefilled Normal Saline Flush Syringe (10ml) contains 0.9% Sodium Chloride USP. Compatible with all Luer-lock needle-free IV access systems. Single-use prefilled syringes help reduce medication errors and contamination risks associated with multi-dose vials. Sterile and for single use only.

Item # PF1010
Packing: 30 Prefilled Syringes / box
Price: \$25.00 /box

Sterile Disposable Syringes



KBMED Sterile Disposable Syringes (without needle) are single-use, individually wrapped syringes designed for precision and ease of use. The ultra-clear barrel features well-defined markings to ensure accurate dosing. These luer lock syringes have a smooth surface for effortless dispensing and an extended gasket to minimize dead space.

Size:	Item#	Packing:	Price:
1ml	#DS3P-1	100 pcs	\$18.00
3ml	#DS3P-3	100 pcs	\$18.00
5ml	#DS3P-5	100 pcs	\$18.00
10ml	#DS3P-10	100 pcs	\$18.00
30ml	#DS3P-30	50 pcs	\$18.00
60ml	#DS3P-60	25 pcs	\$18.00

Micro Syringes with Ultra-Fine Needle



MERITPOINT Micro Syringes with Ultra-Fine Needle provide superior injection comfort and precision. Ideal for injecting medications such as insulin, Botox, collagen, hyaluronic acid, and mesotherapy solutions.

Size:	Item#
0.3ml (1/3cc) 31G x 8mm	#SN03-3108
1.0ml (1cc) 31G x 8mm	#SN10-3108

Packing: Box of 100
Price: \$18.00 /box

Hypodermic Disposable Needles



EXEL® Hypodermic Disposable Needles are engineered for precision and dependability. Made from high quality 304 stainless steel, they deliver a sharp, smooth injection experience. ISO color-coded luer lock hub provides easy identification and a secure fit. Individually sealed and sterilized.

Size:	Item#
18G x 1.5" Blunt Tip	#N-305180
18G x 1" (25mm)	#N-26419
19G x 1" (25mm)	#N-26424
19G x 1.5" (38mm)	#N-26438
20G x 1" (25mm)	#N-26414
22G x 1" (25mm)	#N-26411
23G x 1" (25mm)	#N-26408
30G x 0.5" (13mm)	#N-31000
30G x 1" (25mm)	#N-26439

Packing: Box of 100
Price: **\$15.00 /box**


Blood Collection Tube Holders



METIER Blood Collection Tube Holders are made of a clear plastic for easy viewing of blood collection while providing a secure connection to blood collection sets, needles, and tubes when performing a direct draw into a blood collection tube. Designed to improve safety and compliance during routine blood collection procedures. Non-sterile and for single use only.

Item # **BCH-250**
Packing: 250 pcs /bag
Price: **\$25.00 /bag**

Safety Hypodermic Needles



KBMED Safety Hypodermic Needles are designed for compatibility with all syringes featuring a luer lock or luer slip tip, ensuring seamless clinical injections. Engineered for of use, they offer clear visibility of the injecti site and help prevent accidental needle stick injuries. With a simple one-handed activation, the needle is securely covered immediately after use, providing instant protection. Each needle is sterilized using EO gas and for single use only.

Size:	Item#
21G x 1.5" (38mm)	#SN1B2138
22G x 1" (25mm)	#SN1B2225
23G x 1" (25mm)	#SN1B2325
25G x 1" (25mm)	#SN1B2525

Packing: 100 pcs /box
Price: **18.00 /box**

Luer-to-Luer Fluid Dispensing Connectors



B. BRAUN® Luer-to-Luer Fluid Dispensing Connector is designed for safe and efficient fluid transfer. It enables secure connections between Luer and Luer lock fittings to facilitate liquid transfer or mixing of injectable solutions between syringes while maintaining sterility.

Item # **FDC1000**
Packing: 100 pcs /box
Price: **\$180.00 /box**

Ultra-thin Wall Nano Cannulas



Ultra-thin Wall Nano Cannulas by JBP are designed to reduce patient discomfort, minimize hematomas and limit swelling. Flexible cannula with blunt tip slides smoothly under the skin without causing trauma. Ideal for mesotherapy, placenta injection, Botox injection, collagen injection and local anesthetics. Each cannula comes packaged with a sharp punching needle. Sterilized by EO Gas.

Size:	Item#
22G x 2" (50mm)	#CANN-2250
22G x 2 3/4" (70mm)	#CANN-2270
25G x 1 1/2" (40mm)	#CANN-2540
25G x 2" (50mm)	#CANN-2550
25G x 2 3/4" (70mm)	#CANN-2570
27G x 1 1/2" (40mm)	#CANN-2740
27G x 2" (50mm)	#CANN-2750
30G x 1" (25mm)	#CANN-3025

Packing: 24 pcs /box
Price: **\$99.50 /box**

Sterile Luer Lock Syringe Tip Caps



MEDIGRATIVE® Sterile Syringe Tip Cap has specialized grooves that securely lock it in place, preserving the sterility of all standard luer lock syringes. Easy to identify white color. Individually sealed and sterilized.

Item # **CA001**
Packing: 200 pcs /box
Price: **\$28.00 /box**

Ultra-thin Wall Nanoneedles

Ultra-thin Wall Nanoneedles by JBP are designed to reduce patient discomfort, minimize hematomas and limited swelling. Using a patented process called Micron Sharpness, the nanoneedles allow injection of higher viscosity solutions through a smaller entry to reduce pain and bleeding. Goal is to minimize the surface area so injectors have minimal impact on soft tissue and more satisfied patients. Sterilized by EO gas.



Size: 30G (4mm)
Item # **N-UTW3004**
Packing: 100 pcs /box
Price: **\$65.00 /box**



Size: 33G (13mm)
Item # **N-UTW3313**
Packing: 100 pcs /box
Price: **\$75.00 /box**



Size: 34G (4mm)
Item # **N-UTW3404**
Packing: 100 pcs /box
Price: **\$85.00 /box**

Size: 34G (8mm)
Item # **N-UTW3408**
Packing: 100 pcs /box
Price: **\$85.00 /box**

Microsurgical Blades

SHARD® Microsurgical Blades are developed for minimally invasive procedures and microsurgery. These finer blades are suitable for use in many surgical disciplines. Made of high-grade carbon or stainless steel with an ultra-sharp cutting edge for clean incisions. High surface finishing ensures smooth texture. Each single-use blade is individually wrapped and sterilized by gamma radiation.

	#A600-M61 12 blades / box Price: \$19.50 /box
	#A600-M62 12 blades / box Price: \$19.50 /box
	#A600-M64 12 blades / box Price: \$19.50 /box
	#A600-M65 12 blades / box Price: \$19.50 /box
	#A600-M67 12 blades / box Price: \$19.50 /box
	#A600-M69 12 blades / box Price: \$19.50 /box

	#A600-M90 25 blades / box Price: \$29.50 /box <i>Ultra-Fine</i>
	#A600-M91 25 blades / box Price: \$29.50 /box <i>Ultra-Fine</i>

Microsurgical Blade Handle



SHARD® Microsurgical Blade Handles feature a 4.5" (115mm) long round handle with serrated thumbs for better control. Made of surgical-grade stainless steel that is fully autoclavable. Compatible with all SHARD microsurgical blades.

#A600-MBH
1 Piece
Price: **\$29.00 /piece**

Premium Surgical Blades

SHARD® Premium+ Surgical Blades are available in 5 popular shapes and are used throughout all surgical disciplines including dermatology, oral surgery, orthopedics and reconstructive surgery. Computer-controlled grinding ensures a consistently sharp edge. Sterilized by gamma radiation.

#10 	ITEM# A600-P10-0
#11 	ITEM# A600-P11-0
#12 	ITEM# A600-P12-0
#15 	ITEM# A600-P15-0
#15C 	ITEM# A600-P15-C

Packing: 100 individually sealed blades per box
Price: **\$34.00 /box**

Suture Removal Scissors



AD SURGICAL Suture Removal Scissors are ideal for suture cutting and removal. The curved, blunt-tipped blades help locate and cut buried sutures without traumatizing adjacent sensitive tissue. Each scissors is 4 3/4" long and is fully autoclavable with a lifetime no-rust guarantee. Made in USA.

Item # **SCR-B115**
Packing: 1 pair
Price: **\$45.00 each**

No. 3 Surgical Scalpel Handles

Swann-Morton No. 3 Surgical Scalpel Handles are designed to accommodate blade sizes 6 through 16, offering full compatibility with widely used surgical and dermaplaning blades including 10, 10A, 10R, 11, 12, 14, 15, and 15A. Crafted from solid stainless steel, this handle is built for durability and professional use, offering easy sterilization and long-term performance in clinical and surgical settings. Made in England.



#3 Fitment Graduated Handle (Short, 4.5")
#CS-073S
 Packing: 1 PC
 Price: **\$18.00**



#3 Fitment Graduated Handle (Long, 8")
#CS-073L
 Packing: 1 PC
 Price: **\$28.00**



#3 Fitment Handle (Slim & Short, 4.5")
#CS-07SM9
 Packing: 1 PC
 Price: **\$18.00 /box**



#3 Fitment Handle (Slim & Long, 6.5")
#CS-077S
 Packing: 1 PC
 Price: **\$28.00 /box**



#3 Fitment Handle (Ergonomic w/ Stainless Steel, 5")
#CS-07SM5B-S
 Packing: 1 PC
 Price: **\$84.00 /box**



#3 Fitment Handle (Ergonomic w/ Plastic, 5")
#CS-07SM5B-E
 Packing: 1 PC
 Price: **\$22.00 /box**

Disposable Biopsy Punch



SHARD® Premium+ Disposable Biopsy Punch is ergonomically designed for routine dermal procedures. The razor-sharp seamless, stainless steel blade ensures safe, clean and easy sample extraction. Features a unique, easy-grip plastic handle with pressure assistance platform for superb control and comfort. Hollow handle allows easy specimen retrieval. Each punch is individually sealed and sterilized by gamma radiation.

- | | |
|------------|-----------------|
| Size 1.5mm | Item #A750-BP15 |
| Size 2.0mm | Item #A750-BP20 |
| Size 2.5mm | Item #A750-BP25 |
| Size 3.0mm | Item #A750-BP30 |
| Size 3.5mm | Item #A750-BP35 |
| Size 4mm | Item #A750-BP40 |
| Size 5mm | Item #A750-BP50 |
| Size 6mm | Item #A750-BP60 |
| Size 8mm | Item #A750-BP80 |

Packing: 10 pcs/box
 Price: **\$26.00 /box**

Disposable Surgical Scalpels

Our economical SHARD® Disposable Scalpels incorporate a razor-sharp SS blade on an ergonomic plastic handle for excellent performance and control. Each scalpel is sealed individually and sterilized by gamma radiation.



Packing: 10 individually sealed scalpels per box
 Price: **\$6.50 /box**

Disposable Dermal Currettes



SHARD® Premium+ Disposable Dermal Currettes are used for scraping skin lesions and growths such as warts or melanoma. It is a single-ended instrument available with a round loop tip. The stainless steel loop tip is sharp on the inside and blunt on the outside so the lesion can be removed, but surrounding skin is not damaged. Each curette is individually sealed and sterilized by gamma radiation.

- | | |
|----------|----------------|
| Size 2mm | Item #A700-DC2 |
| Size 3mm | Item #A700-DC3 |
| Size 4mm | Item #A700-DC4 |
| Size 5mm | Item #A700-DC5 |
| Size 7mm | Item #A700-DC7 |

Packing: 10 pcs/box
 Price: **\$29.50 /box**

Disposable Dermal Biopsy Blade



SHARD® Premium+ Dermal Biopsy Blades are sterile, flexible scalpels designed for cutaneous surgery, including shave biopsies, saucerization of flat lesions, and leveling of pedunculated lesions. Our unique dermal blade design features a comfortable and protective finger guard along with an ultra-sharp, flexible blade. This flexibility allows the blade to be positioned at the correct angle for the intended procedure.

Item # **A750-DBB**
 Packing: 50 pcs/box
 Price: **\$105.00 /box**

SWANN-MORTON Single Use Surgical Blade Remover



The SWANN-MORTON Single Use Blade Remover is engineered to extract surgical blades safely and efficiently from No.3 and No.4 handles, supporting blade sizes across a wide range. Designed for single-use applications, this tool ensures safe handling and disposal during surgical, veterinary, or aesthetic procedures.

Item # **CS-48SM**
 Packing: 100 pcs/box
 Price: **\$32.00 /box**

SWANN-MORTON Surgical Blade Remover Unit



The Swann-Morton Blade Remover Unit is designed for safe and efficient hands-free removal of used surgical blades, making it an ideal solution for clinics, aesthetic practices, veterinary offices, and educational labs. Its simple press-and-eject mechanism ensures quick blade removal without the need to touch the blade directly. Each non-sterile blade remover unit securely holds up to 100 used blades

Item # **CS-48**
 Packing: 10 pcs/case
 Price: **\$28.00 /case**

QLICKSMART BladeFlask EVO Blade Remover



The QlickSmart BladeFlask EVO Blade Remover can accommodate a broad range of scalpel handles and blades, including large autopsy scalpels, thick and round handles, even hexagonal shaped versions – giving you greater freedom in using any scalpel you want. It is simple to use with just one hand. The built-in blade counter indicates when the unit is full. Embrace safety with the new standard in scalpel blade removals.

Item # **CS-48EVO**
 Packing: 1 piece
 Price: **\$50.00**

BladeFlask Mounting Bracket
 Item # **CS-48QSB**
 Packing: 1 piece
 Price: **\$10.00**

Universal Surgical Drapes



AD SURGICAL Sterile Universal Surgical Drapes are made of an absorbent and completely impervious PE/ PP material to prevent strike-through and subsequent contamination of the surgical site. Each drape has a peel-able adhesive tape and can be used alone or in combination, giving flexibility to individual procedural needs. Individually wrapped and sterilized by EO Gas. Latex free.

HALF DRAPE WITH TAPE

Sterile Half Drape with Adhesive Tape

- 60in x 44in (152cm x 112cm)
- Absorbent, Impervious PE/PP Material

Item # **A400-HFD**
 Packing: 10 individually sealed packs per box
 Price: **\$38.00 /box**

MINOR PROCEDURE DRAPE with adhesive aperture

Sterile Drape with Adhesive Tape

- 22in x 25in (56cm x 64cm)
- Aperture Size: 2.5in (6cm) Diameter
- Absorbent, Impervious PE/PP Material

Item # **A400-APD**
 Packing: 20 individually sealed packs per box
 Price: **\$32.50/box**


UTILITY DRAPE WITH TAPE

Sterile Drape with Adhesive Tape

- 22in x 25in (56cm x 64cm)
- Absorbent, Impervious PE/PP Material

Item # **A400-UTD**
 Packing: 20 individually sealed packs per box
 Price: **\$32.50/box**

Sterile Hose & Cable Covers



EZ DRAPE® Sterile Hose and Cable Covers are user friendly and effective. Simply insert the handpiece or other surgical instrument through the ring and pull the plastic sleeve down. Sterile hose cover helps maintain a clean and safe surgical environment. Individually wrapped and sterilized by EO Gas. Latex free.

Ring diameter: 2.2" Length: 60" long

Item # **A400-1000**
 Packing: 24 individually sealed packs per box
 Price: **\$64.00/box**

Latex Surgical Gloves (Powder-Free)



SUPERFIT® Sterile Latex Surgical Gloves are made from ergonomic molds that help reduce hand fatigue during operation. Powder-free reduces the risk of latex allergy. Straight finger design and wet donning ability make it very easy to don either with dry or damp hands. Each pair of Left/ Right fitted gloves are sealed in a sterile pouch and sterilized by gamma / ebeam radiation.

Sizes	Item#
Size 6.0	A520-LTX60
Size 6.5	A520-LTX65
Size 7.0	A520-LTX70
Size 7.5	A520-LTX75
Size 8.0	A520-LTX80

Packing: 50 pairs per box
 Price: **\$54.00 /box**

Absorbable Hemostatic Dressings Silicone Scar Removal Sheets



BENACEL® Absorbable Hemostatic Dressing is a specially formulated gauze made of biocompatible oxidized cellulose and contains no harmful animal by-products or chemical additives. When applied to a wound, it produces hemostasis almost instantly by forming a gel-like substance upon contact with blood. BENACEL® is used to control bleeding during surgical procedures, or as a topical dressing on superficial open wounds. Made in USA.

0.5cm x 0.7cm plug
Item # **B01-0108**

Packing: 8 pcs /box
Price: **\$42.00 /box**

0.6cm x 0.8cm plug
Item # **B01-0308**

Packing: 8 pcs /box
Price: **\$46.00 /box**

0.6" x 0.6" (1.5cm x 1.5cm)
Item # **B01-0212**

Packing: 12 pcs /box
Price: **\$55.00 /box**

2" x 2" (5cm x 5cm)
Item # **B01-0510**

Packing: 10 pcs /box
Price: **\$66.00 /box**

2" x 4" (5cm x 10cm)
Item # **B02-1010**

Packing: 10 pcs /box
Price: **\$148.00 /box**

Collagen Wound Dressing Pads



Made of 100% bovine collagen, WOUND FREE® Resorbable Collagen Wound Dressing Pads encourage the body to heal wounds quickly and effectively. Once applied to the wound, the collagen dressing forms a soft and comfortable gel on contact with the fluid that promotes healing throughout the day. For medical and dental use. Sterile and for single use only.

2" x 2" (5cm x 5cm)
Item # **CWD-0505**

Packing: 10 pcs /box
Price: **\$75.00 /box**



WOUND FREE® Silicone Scar Removal Sheets contain medical-grade silicone that helps hydrate scar tissue and mimics the natural barrier function of healthy skin to make scar softer, flatter and less visible. Each self-adherent sheet can last for up to 7-10 days of use. Color: Tan.

#SCS-T1R 1.5" x 2.4" (4cm x 6cm) 10 PCS / Box

Price: **\$16.00 / box**

#SCS-T5R 1.2" x 6" (3cm x 15cm) 10 PCS / Box

Price: **\$29.50 / box**

Silicone Foam Dressings



WOUND FREE® Sterile Silicone Foam Dressings combine the superior absorption of foam and the gentle adhesion of silicone. The foam pad absorbs large amount of exudates. The silicone adhesive border provides gentle fixation, allowing the dressing to be removed or reapply when needed. Suitable for mild to moderately exuding wounds.

FMD-S4 4"x4" (10cm x 10cm) 10PCS/Box

Price: **\$22.00 / box**

FMD-S6 6"x6" (15cm x 15cm) 10PCS/Box

Price: **\$36.00 / box**

Adhesive Island Dressings



WOUND FREE® Sterile Adhesive Island Dressing is a gentle soft cloth bordered gauze dressing indicated for wounds with light to moderate exudation. One-step application is an excellent alternative to tape and gauze to help save user time. The non-adherent central gauze pad provides absorption of light to moderate drainage.

ISD-P0708 2.75"x 3.15" (7cm x 8cm) 55PCS/Box

Price: **\$12.00 / box**

ISD-P1010 4"x 4" (10cm x 10cm) 40PCS/Box

Price: **\$14.00 / box**

Super Absorbent Dressings



WOUND FREE® Sterile Super Absorbent Dressing is a new generation of advanced wound dressing. It can be used as either a primary or secondary dressing with the purpose of using its extremely absorbent core for exudate management. The self-adherent super absorbent polymer dressing core effectively soaks up high amounts of wound exudate from pressure and venous ulcers, surgical incisions, burns, donor grafts and chronic wounds.

Self-Adherent Version

SAD-AS4 4"x4" (10cm x 10cm) 10PCS/Box

Price: **\$12.00 / box**

Reusable Cold Gel Packs



AD SURGICAL Reusable Cold Gel Pack helps relieve pain and swelling after a surgical procedure or during the initial stage of an acute injury. 4" diameter with soft cloth backing provides comfortable cold therapy for up to 20 minutes. Cold gel pack is commonly used for dermal fillers, chemical peels, facial injections, oral surgery and other surgical procedures.

#T430-BU 4" Round 10 PCS / Bag

Price: **\$25.00 / Bag**

Soft Stretch Facial Wraps



AD SURGICAL Post-Op Soft Stretch Facial Wraps by COOL JAW® feature an exclusive, hand-free patented design. Wraps are made with a spandex fabric that is extremely soft on the patient's skin. They feature a continuous outer sleeve that extends cheekbone to cheekbone to ensure full coverage and compression over the surgical area. This sleeve is specifically designed to fit our gel packs and allows patients to easily insert (and replace) gel packs without having to remove the wrap. The wrap features two rows of Velcro closure at the top of the head to provide a universal, adjustable fit.

Face Wrap With Chin Cup, Black

Item# **W800C-BLK**

Price: **\$16.00 /piece**

Face Wrap Without Chin Cup, White

Item# **W800-WHT**

Price: **\$16.00 /piece**

Post-op Therapy Gel Packs

AD SURGICAL Hot/Cold Therapy Gel Packs are frequently used by cosmetic and plastic surgery patients to reduce swelling and bruising. They are also used by patients after oral surgeries including wisdom teeth extractions, reconstructive jaw surgery, and implant surgery. Our reusable Hot/Cold Gel Packs stay pliable when frozen and remain cold for over an hour.



Cold Gel Pack
3" x 5" Clear

Item# **G400-CLR**

10 PCS
Price: **\$23.00**



Hot + Cold Gel Pack
3" x 6" Galaxy Blue

Item# **GB425-BL**

12 PCS
Price: **\$32.00**

Sterile Oral Surgery Drape Packs



O-PAK® Sterile Oral Surgery Drape Pack helps maintain a clean and safe surgical environment to prevent postoperative infections and cross contamination. Ideal for dental implant surgery or other aseptic procedures. Simplify procedural setup with 14 essential disposable products included in one convenient package. Sterilized by EO Gas. Latex free.

Each pack contains all the items listed below :

(Items are not sold separately)

1. 2 Procedure gowns, universal size
2. 2 Facemasks with ties, 2 Bouffant caps
3. Absorbent cotton scrub towels, 20 x 40 cm, 1 Surgical instrument tray
4. 1 Impervious patient drape with tape, 112 x 152 cm
5. 2 Table or mayo stand covers with adhesive tape, 64 x 43 cm
6. 1 Absorbent O/R towel
7. 1 Saliva ejector, 1 Surgical aspirator tip, 2 Cotton tip applicators, 25 Gauze swabs
8. 4 Light handle/ chair button barrier films
9. 1 Medical wrapping sheet, 90 x 90 cm

Item # **A200-3000**

Packing: 6 individually sealed packs per box

Price: **\$128.00 /box**

Oral Tissue Adhesives

PERIACRYL® Oral Tissue Adhesive is a biocompatible and non-toxic cyanoacrylate tissue adhesive that is useful in securing sutures, stabilizing membranes, covering bone grafting material post-extraction, or dressing donor and recipient sites. Violet coloring allows for better visibility. Available in multi-use and single-use versions.

Multi-Use Kit (High Viscosity)

Item# **GLU-PERI50**



Packing: 5 ml bottle, 50 pipettes + 1 Tray

Price: **\$130.00 /box**

Single-Use 12 Pack (High Viscosity)

Item# **GLU-PERI12**



Packing: 12 x (0.2 ml) / box

Price: **\$130.00 /box**

Universal Implant Drill Stopper Set Implant Impression Trays



Autoclavable

EZ STOPPER® Universal Implant Drill Stopper Set is a set of medical grade, autoclavable silicone tubings. Each tubing can be accurately measured and cut to the desired length in seconds using the built-in lockable metal stop in the autoclavable storage tray. The trimmed silicone stoppers will stretch-fit over most surgical drills to provide precise drill depth control during dental implant placement. Each silicone stopper is color-coded to fit the corresponding surgical drill sizes as follows:

- RED Ø 1.8 – 2.4mm
- GREEN Ø 2.5 – 3.4mm
- YELLOW Ø 3.5 – 4.4mm
- BLUE Ø 4.5 – 5.6mm

EZ Stopper Starter Kit:


Item #D017-STPS
Packing: 1 set (includes 4 color coded tube)
Price: \$160.00 /set

Silicone Tubing Refills:




- Red** (Item #D017-STP18)
- Green** (Item #D017-STP25)
- Yellow** (Item #D017-STP35)
- Blue** (Item #D017-STP45)

Packing: 4 PCS per refill pack
Price: \$28.00 /pack



IMPLANT TRAY® Implant Impression Tray has an innovative patented design that features a wide array of clear and easily removable panels to accommodate open-tray impression techniques. This eliminates guesswork and minimizes time needed to drill out holes in the impression tray. Made in USA.



- Large Upper** Item# D020-U3L
- Medium Upper** Item# D020-U2M
- Small Upper** Item# D020-U1S
- Large Lower** Item# D020-L3L
- Medium Lower** Item# D020-L2M
- Small Lower** Item# D020-L1S

Packing: 12 trays per box
Price: \$30.00 /box

Assortment Starter Package

Item # D020-COMBO
Contents: 6 of each sizes (Large Upper, Medium Upper, Small Upper, Large Lower, Medium Lower, Small Lower)

Packing: 36 trays per box
Price: \$84.00 /box

Implant Abutment PTFE Fillers



Autoclavable

i-PLUG® Implant Abutment PTFE Filler is a packable filler material for sealing abutments. The one-size PTFE rod (Ø1.8mm x 50mm) can be cut to custom fit the abutment hole of any screw-retained or cement-retained crown. i-PLUG® is fully autoclavable to prevent cross-contamination.

Item # D030-050
Packing: 1 box (15 pcs)
Price: \$49.50 /box

Filling & Removal Instrument



Double-sided stainless steel instrument with both a plugger end and a barbed retrieval end. Designed for easy placement or removal of our i-PLUG PTFE abutment fillers. Fully autoclavable & rust resistance.

Item # D030-017
Packing: 1 double-sided instrument
Price: \$39.50

AUGMA Degranulation Kit



AUGMA Degranulation Kit is intended for removing granulation tissue before augmentation. Ideal for removing all remnants of soft tissue which are well attached to the bone and grab them out without harming the bone. Each kit includes the following bur sizes: 1mm, 2.5mm, 3mm, 3.5mm.

Item # D035-D-KIT
Packing: 1 Kit
Price: \$230.00 /Kit

Thermoplastic Surgical Guides for Pilot-Drilling

EZ STENT® Thermoplastic Surgical Guide for Pilot-Drilling is our original patented design that allows fabrication of an accurate surgical guide in minutes using hot water alone. The template has a titanium drill sleeve with inner diameter 2.3mm that is compatible with the pilot drill of most dental implant systems. EZ STENT® should be cold sterilized prior to usage and is for single use only.

SINGLE TOOTH SET

EZ STENT® Single Tooth Set allows accurate preparation of an osseous site in a tooth-bound region. Titanium drill sleeve with inner Ø2.3mm is compatible with the pilot drill of most dental implant systems.



Item # D015-003-P
Contents:
 5 Thermoplastic templates for Single-Tooth sites
 - Sleeve height 10mm, sleeve inner Ø2.3mm
 5 Guide pins (stainless steel), Length 20mm, Ø2.3mm
 1 Drill for dental laboratory (steel), Ø2.3mm
Price: \$160.00 /set

FREE END SET

EZ STENT® Free End Set allows accurate preparation of an osseous site in the most distal region. Titanium drill sleeve with inner Ø2.3mm is compatible with the pilot drill of most dental implant systems.



Item # D015-004-P
Contents:
 5 Thermoplastic templates for Free-End Situations
 - Sleeve height 10mm, sleeve inner Ø2.3mm
 5 Guide pins (stainless steel), Length 20mm, Ø2.3mm
 1 Drill for dental laboratory (steel), Ø2.3mm
Price: \$160.00 /set

2-UNIT BRIDGE SET

EZ STENT® 2-Unit Bridge Set allows the surgeon to accurately locate and drill the holes of two consecutive implants with optimal spacing, angulation and parallelism. Titanium drill sleeves with inner Ø2.3mm are compatible with the pilot drill of most dental implant systems.



Item # D015-008-2U
Contents:
 4 Thermoplastic templates (titanium/ polymer)
 - 2 Tooth-Bound and 2 Free End Templates
 - Sleeve height 6mm, sleeve inner Ø2.3mm
 8 Guide pins (stainless steel), Length 20mm, Ø2.3mm
 1 Drill for dental laboratory (steel), Ø2.3mm
Price: \$160.00 /set

Pilot Surgical Drills (for EZ Stent®)

Our Ø2.3mm Pilot Surgical Drills are perfect fit for the EZ STENT® Pilot-Drill Sets. Feature wide-band, easy to see, laser-etched depth markings at 8mm, 10mm, 11.5mm, 13mm and 16mm. Made in USA.



Item # D015-023L
Content:
 1 Long pilot surgical drill (stainless steel)
 Length: 38mm
Price: \$68.00 /drill

Item # D015-023S
Content:
 1 Short pilot surgical drill (stainless steel)
 Length: 30mm
Price: \$68.00 /drill

Synthetic Bioactive Bone Grafts



UNIGRAFT® Bioactive Bone Graft is made of a synthetic bioactive glass material. It offers an osteogenic framework for new bone growth and bolsters a higher level of osteoblastic activity to enhance bone regeneration. UNIGRAFT® is completely resorbable and transforms into natural bone with time. Made in USA.

Bone Graft (200-600um)

Packing: 5 vials per box

0.5 gram Item# U03-0105
Price: **\$182.00 /Box**

1.0 gram Item# U03-0205
Price: **\$292.00 /Box**

2.0 gram Item# U03-0305
Price: **\$396.00 /Box**

BOND APATITE Bone Graft Cements



BOND APATITE by Augma Materials is a composite graft, made of Biphasic Calcium Sulfate mixed with HA granules in a controlled particle size distribution, intended to fill or augment a broad range of bony defects. It fully resorbs and is replaced by bone during the healing process. It allows minimal invasive surgical procedure, reduced treatment time and easy manipulation of the graft material for the clinician.

Item# HA-BB-PACK

Packing: 1cc Syringe / Box

Price: **\$135.00 /Box**

Mineralized Allograft In Syringes



OSSEOSEAL™ bone in syringe is a mixture of allograft mineralized cortical and cancellous powder in a convenient syringe with a filter cap! Curved syringe for better access to site. The cancellous and cortical blend forms an osteoconductive scaffold providing volume enhancement and effective site development for successful ridge preservation and dental implant placement. Sterile and for single use only.

Bone Graft (250-800um)

0.3cc Item# OS-03CC
Price: **\$44.00 /Pack**

0.5cc Item# OS-05CC
Price: **\$66.00 /Pack**

1.0cc Item# OS-10CC
(2 of 0.5cc syringes in 1 Pack)
Price: **\$100.00 /Pack**

Attachable Intraoral Wound Dressings



AUGMA SHIELD® is an adhesive intraoral patch used to protect a wound such as an extraction socket, tissue grafting & implant site, and ulcers. Easy to use, just peel off and stick on the wound site. The non-resorbable polymer helps cover the affected area until soft tissue closes and healing begins.

25mm x 15mm Intraoral Dressing

Item# AS-20PER

Packing: 20 Pieces / Box

Price: **\$99.00 /Box**

Bovine Bone Grafts



OSSEOSEAL™ Bovine Bone Graft is a xenograft produced by removal of organic components from bovine bone. It provides a supportive structure for osteoconduction. The porous bone mineral matrix has macro and microscopic structures that mimic human bone. The formation and ingrowth of new bone at implantation site is favored, due to its trabecular architecture, interconnecting macro and micro pores.

Bone Graft (200-1000um)

Packing: 1 vial per box

0.5 gram Item# OSX-05CC
Price: **\$38.50 /Box**

1.0 gram Item# OSX-10CC
Price: **\$55.00 /Box**

4.0 gram Item# OSX-40CC
Price: **\$132.00 /Box**

Allograft Bone Fiber Plugs



REVOSS® Matrix Bone Fiber Plugs are composed of a human allograft mixture of demineralized cortical bone fiber and mineralized cancellous bone particulates. The unique, compressed cylinder shape of REVOSS Matrix allows for easy insertion and expansion within the bone void, eliminating the need for particulate mixing or membranes to facilitate bone regeneration.

11mm x 7.5mm Bone Insert

Item# RM1-0101

Packing: 1 Piece / Box

Price: **\$96.00 /Box**

Mineralized Allograft in Vials



REVOSS®/ OSSEOSEAL® Mineralized Human Allograft in Vials is a sterile processed human allograft tissue ideal for a variety of surgical applications. It comprises 50/50 mineralized cortical / cancellous bone powder. All tissue was recovered, processed, stored and distributed for use in accordance with the standards of AATB and FDA requirements.

- 0.25 CC Item# RB5-1025 Price: **\$46.00 /Box**
- 0.5 CC Item# RB5-1050 Price: **\$75.00 /Box**
- 1.0 CC Item# RB5-1100 Price: **\$90.00 /Box**
- 2.5 CC Item# OEU0250 Price: **\$135.00 /Box**
- 5.0 CC Item# OEU0250-2X Price: **\$270.00 /Box**

Demineralized Allograft Putty



REVOSS™ Demineralized Allograft Putty is a demineralized bone matrix (DBM) allograft made of donated human tissue. It is packaged in an easy-to-use prefilled syringe for direct application with no need for pre-treatment in a variety of surgical applications. All tissue was recovered, processed, stored and distributed for use in accordance with the standards of AATB and FDA requirements.

Demineralized Allograft Putty in Syringes

- 0.7 cc Item# RP2-0070 Price: **\$136.00 /1pc**
- 1.0 cc Item# RP2-0100 Price: **\$158.00 /1pc**
- 2.5 cc Item# RP2-0250 Price: **\$370.00 /1pc**

Resorbable Bone Substitutes



BICERA® Resorbable Bone Substitute is a bioceramic material that is composed of 60% hydroxyapatite (HAP) and 40% beta-tricalcium (β-TCP). It can be used in common dental procedures such as sinus lifting, periodontal bone fillings, and filling of peri-implant defects. BICERA provides a bone void filler that resorbs and is replaced by bone during the healing process.

Bicera HA/TCP Granule 1-2mm

- 2cc Item# BC-G22 Price: **\$160.00 /Box**
- 5cc Item# BC-G02 Price: **\$240.00 /box**
- 10cc Item# BC-G07 Price: **\$400.00 /Box**



Bicera Block
(5mm x 5mm x 20mm)
Item# BC-B01
Packing: 4 Pcs
Price: **\$220.00 /Box**

Resorbable Collagen Plugs



HEALIAID® Resorbable Collagen Plugs are designed to deliver dependable wound dressing at extraction and biopsy sites. Shaped for ease of use and compressible to suit open wounds, these resorbable dressings are inexpensive and easy to work with. When used to dress fresh extraction sockets, the collagen plugs protect exposed surfaces from external irritants and further injury. The Type-1 purified bovine collagen is non-pyrogenic, and typically reabsorbed by the body within 2 to 4 weeks after placement. Each collagen plug comes sterile and is ready for use without preparation.

10mm x 20mm Collagen Plug

Item# HAD001020

Packing: 10 Plugs / Box

Price: **\$105.00 /Box**

Bilayer Collagen Matrix



HEALIAID® Bilayer Collagen Matrix is designed for soft tissue regeneration and sealing extraction sockets in the oral cavity. It is composed of two layers: the dense outer layer provides stability while allowing open healing; the spongy inner layer supports blood clot stabilization and promotes early vascularization and good soft-tissue ingrowth.

10mm x 10mm

(Small- for sealing extraction sockets)

Item# HADB101030

Packing: 5 Pieces / Box

Price: **\$220.00 /Box**

15mm x 15mm

(Medium- for sealing molar extraction sockets)

Item# HADB151530

Packing: 5 Pieces / Box

Price: **\$280.00 /Box**

25mm x 30mm

(Large- for soft tissue regeneration)

Item# HADB253030

Packing: 1 Piece / Box

Price: **\$140.00 /Box**

30mm x 40mm


(X-Large- for soft tissue regeneration)

Item# HADB304030

Packing: 1 Piece / Box

Price: **\$200.00 /Box**

Resorbable PLA / PGA Barrier Membranes



CYTOFLEX RESORB® Resorbable PLA / PGA Barrier Membrane consists of three integral layers designed to promote gingival tissue attachment and prevent fibroblast down growth. Made of biocompatible poly (lactide-co-glycolide) copolymers, CYTOFLEX RESORB® is non-pyrogenic, non-immunogenic and maintains a barrier framework for over two months after implantation. The resorbable membrane does not require pre-soaking before administration, and completely dissolves within six months. Made in USA.

12mm x 24mm Resorbable Membrane

Item# C03-0301

Packing: 1 per pack

Price: **\$60.00 /pack**

20mm x 25mm Resorbable Membrane

Item# C03-0201

Packing: 1 per pack

Price: **\$88.00 /pack**

30mm x 40mm Resorbable Membrane

Item# C03-0101

Packing: 1 per pack

Price: **\$124.00 /pack**

ePTFE Barrier Membranes



CYTOFLEX TEF-GUARD® ePTFE Barrier Membrane is a non-resorbable membrane made of microporous ePTFE material. It is designed to resist the penetration of fibroblasts and bacteria, while simultaneously allowing the exchange of interstitial fluids through the membrane. This enhances gingival tissue attachment and provides a favorable protected environment for neo-vascularization and re-population of osteoblasts in the bony defects. Made in USA.

12mm x 24mm ePTFE Membrane

Item# C04-0101

Packing: 1 per pack

Price: **\$4.00 /pack**

25mm x 30mm ePTFE Membrane

Item# C04-0201

Packing: 1 per pack

Price: **\$72.00 /pack**

Resorbable Collagen Membranes



COLLAFLEX™ Resorbable Collagen Membrane is derived from purified porcine pericardium. The non-crosslinked membrane retains the natural microfibril structure of the pericardium, providing excellent biocompatibility and extended barrier function for guided bone regeneration applications. Highly drapable and can be easily fixated with sutures and bone tacks.

15mm x 20mm Collagen Membrane

Item# C07-0101

Packing: 1 per pack

Price: **\$75.00 /pack**

20mm x 30mm Collagen Membrane

Item# C07-0201

Packing: 1 per pack

Price: **\$100.00 /pack**

30mm x 40mm Collagen Membrane

Item# C07-0301

Packing: 1 per pack

Price: **\$166.00 /pack**

Titanium-Reinforced ePTFE Barrier Membranes

CYTOFLEX Ti-ENFORCED® Titanium-reinforced ePTFE Barrier Membrane is a multi-layer, non-resorbable membrane. It is a barrier material with a titanium frame enclosed within two layers of ePTFE material. The titanium reinforcement is intended for space and shape maintenance to contain bone grafting material and to minimize graft movement during wound healing. The Ti-Enforced™ membranes are textured on both surfaces to improve handling and adhesion between tissue and membrane. Made in USA.



12mm x 24mm

Item# C05-0301

Packing: 1 per pack

Price: **\$109.00 /pack**



14mm x 24mm

Item# C05-0501

Packing: 1 per pack

Price: **\$136.00 /pack**



19mm x 26mm

Item# C05-0901

Packing: 1 per pack

Price: **\$148.00 /pack**



20mm x 25mm

Item# C05-1301

Packing: 1 per pack

Price: **\$180.00 /pack**



25mm x 30mm

Item# C05-1701

Packing: 1 per pack

Price: **\$200.00 /pack**



32mm x 40mm

Item# C05-4101

Packing: 1 per pack

Price: **\$245.00 /pack**

Sterile Oral Exam Packs



STERILE i-PAK® is an economical 3-in-1 intraoral exam pack that eliminates cross-contamination, improves efficiency and saves on sterilization costs. Each i-pak® includes 1 disposable dental mirror with anti-fog coating and 1 patented disposable double-ended Explorer & Probe. Probe is marked 3-6-9 mm universal style, tapered with round tip for excellent readability and patient comfort. Safe for use on titanium implants. The explorer is made with surgical grade stainless steel with an ultra-sharp tip. Good for routine exam, emergency treatment, off-site visit and high volume oral exam performed by schools, government and other agencies. Sterilized by EO Gas. Latex free.

Item # D002-009-P

Packing: 50 individually sterilized packs /box

Price: **\$98.50 /box**

Titanium Mesh Membranes

CYTOFLEX® Titanium Mesh Membrane is made of a pure, implantable titanium mesh with precision profile. It is designed for use to ensure three-dimensional reconstruction of alveolar bone defects and to facilitate bone augmentation through adequate fixation of the augmentation material. Made in USA.



12mm x 25mm

Mesh Non-Resorbable Membrane

Item# C02-4101

Packing: 1 per pack

Price: **\$48.00 /pack**



20mm x 25mm

Mesh Non-Resorbable Membrane

Item# C02-4301

Packing: 1 per pack

Price: **\$60.00 /pack**



25mm x 30mm

Mesh Non-Resorbable Membrane

Item# C02-4401

Packing: 1 per pack

Price: **\$60.00 /pack**

AD Surgical
1296 Kifer Road, Suite 608,
Sunnyvale CA 94086 USA

FREE SAMPLES
(for selected products)



AD Surgical

Shop online: www.AD-SURGICAL.com

Order by phone: 888.841.8481

A Simple Promise

Our promise is to deliver the best customer experience by offering:

- Same day shipping on most orders
- A 90-day risk free return policy on all products
- A team of friendly and knowledgeable customer support specialists to ensure your complete satisfaction

We look forward to serving you soon.



Scan to learn more about us

Please keep this catalog for reference.

All prices are in U.S. dollars and are subject to applicable taxes according to state and local laws. Availability, prices and delivery rates are subject to change without notice. Please visit our website for the most updated pricing & product details.