



# AD Surgical

*Top Surgical Essentials For Less Since 2000*

**JAN 2025**



*Expertly serving DENTAL, MEDICAL & VETERINARY PROFESSIONALS for over 25 years*

888.841.8481

[www.AD-Surgical.com](http://www.AD-Surgical.com)

# CONTENTS


01	Sutures	Page 02
02	Skin Marking & Closure	Page 17
03	IV & Blood Collection	Page 19
04	Needles and Sharps	Page 21
05	Wound Care	Page 24
06	Infection Control	Page 26
07	Oral & Implant Surgery	Page 28
08	Dental Implant Accessories	Page 29
09	Regenerative Biomaterials	Page 32
10	Dental Exam & Screening	Page 37

## We offer 4 versions of surgical sutures:



### UNIFY® Ultra-Premium Microsutures

UNIFY® Ultra-Premium Microsutures are made with a special reinforced alloy to resist needle bending and distortion during surgery. Advanced micro taper allows minimal force during tissue penetration. These extra-fine sutures are commonly used in microsurgeries such as nerve repair or microvascular surgeries.



### UNIFY® Premium+ Surgical Sutures

UNIFY® Premium+ Surgical Sutures offer our highest performing suture needles for the most demanding surgical procedures. Precision Point Needles are electropolished for increased sharpness. Our proprietary UNIALLOY™ reinforced stainless steel and MULTICUT™ needle technology helps maintain sharpness and shape up to 3X longer than regular needles.















### UNIFY® Surgical Sutures








UNIFY® Surgical Sutures are top quality sutures at factory-direct prices than are up to 80% lower than our competitors. Made of high grade stainless steel with an ultra-sharp needle point to provide strength, flexibility and atraumatic tissue penetration. With UNIFY® sutures, you'll never choose between price and quality.



### UNIFY® Veterinary Surgical Sutures

UNIFY® Veterinary Surgical Sutures offer premium quality at affordable prices, designed for the most demanding surgical procedures. Crafted from high-grade 300 series stainless steel, these sutures feature ultra-sharp needle points for exceptional strength, flexibility, and smooth, atraumatic tissue penetration.





Needle Size	Needle Shape	Thread Size (USP)	Thread Length	Thread Color	Item #	Qty.	Price/box
<b>NYLON MICROSUTURES</b> <span style="float: right;">Comparable to Ethilon® * and Monosof® **</span>							
<b>SMALL</b> 13mm, 3/8 Circle		7/0	18" (45cm)	Black	S-N718SP13	12 pcs	\$69.50
<b>MICRO</b> BV100-3 (MV-70-3) 4mm, 3/8 Circle		10/0	5" (13cm)	Black	XXS-N1005T4	12 pcs	\$225.00
<b>MICRO</b> BV130-4 (MV-100-4) 5mm, 3/8 Circle		9/0	5" (13cm)	Black	XXS-N905T5	12 pcs	\$225.00
<b>MICRO</b> BV130-5 (MV-135-5) 6.4mm, 3/8 Circle		8/0	7" (18cm)	Black	XXS-N807T6	12 pcs	\$145.00
<b>MICRO</b> TG-140-8 (SE-140-8) 6.4mm, 3/8 Circle		10/0	7" (18cm)	Black	XXS-N1007S6	12 pcs	\$145.00
<b>MICRO</b> TG-140-8 (SE-140-6) 6.4mm, 3/8 Circle		10/0	12" (30cm)	Black	XXS-N1012S6D	12 pcs	\$169.50
<b>PREMIUM NYLON</b> <span style="float: right;">Comparable to Ethilon® * and Monosof® **</span>							
<b>LARGE</b> PS-1 / P-14 / PC33 24mm, 3/8 Circle		2/0	27" (69cm)	Black	PLN-227R24	12 pcs	\$43.50
		3/0	27" (69cm)	Black	PLN-327R24	12 pcs	\$43.50
		4/0	27" (69cm)	Black	PLN-427R24	12 pcs	\$43.50
<b>MEDIUM</b> PS-2 / P-12 / PC31 19mm, 3/8 Circle		3/0	27" (69cm)	Black	PMN-327R19	12 pcs	\$43.50
		4/0	27" (69cm)	Black	PMN-427R19	12 pcs	\$43.50
		3/0	18" (45cm)	Black	PMN-318R19	12 pcs	\$39.00
		4/0	18" (45cm)	Black	PMN-418R19	12 pcs	\$39.00
<b>MEDIUM</b> PS-3 / P-11 / PC34 16mm, 3/8 Circle		4/0	18" (45cm)	Black	PMN-418R16	12 pcs	\$39.00
		5/0	18" (45cm)	Black	PMN-518R16	12 pcs	\$39.00
<b>SMALL</b> P-3 / P-13 / C3 13mm, 3/8 Circle		4/0	18" (45cm)	Black	PSN-418R13	12 pcs	\$39.00
		5/0	18" (45cm)	Black	PSN-518R13	12 pcs	\$39.00
		6/0	18" (45cm)	Black	PSN-618R13	12 pcs	\$39.00
<b>NYLON</b> <span style="float: right;">Comparable to Ethilon® * and Monosof® **</span>							
<b>XX-LARGE</b> KS / SC-2 60mm		2/0	30" (75cm)	Black	XXL-N230K60	12 pcs	\$28.00
<b>XX-LARGE</b> CT / GS-24 / T10 40mm, 1/2 Circle		2	30" (75cm)	Black	XXL-2N30T40	12 pcs	\$28.00

Needle Size	Needle Shape	Thread Size (USP)	Thread Length	Thread Color	Item #	Qty.	Price/box
<b>NYLON</b> <span style="float: right;">Comparable to Ethilon® * and Monosof® **</span>							
<b>X-LARGE</b> CT-1 / GS-21 / T8 36mm, 1/2 Circle		2/0	30" (75cm)	Black	XL-N230T36	12 pcs	\$28.00
<b>X-LARGE</b> FSL / C-16 30mm, 3/8 Circle		0	30" (75cm)	Black	XL-N030R30	12 pcs	\$28.00
<b>LARGE</b> FS-1 / C-14 / C7 24mm, 3/8 Circle		2/0	30" (75cm)	Black	L-N230R24	12 pcs	\$28.00
		3/0	30" (75cm)	Black	L-N330R24	12 pcs	\$28.00
		0	18" (45cm)	Black	L-N018R24	12 pcs	\$28.00
<b>MEDIUM</b> FS-2 / C-13 / C6 19mm, 3/8 Circle		3/0	30" (75cm)	Black	M-N330R19	12 pcs	\$28.00
		4/0	30" (75cm)	Black	M-N430R19	12 pcs	\$28.00
		3/0	18" (45cm)	Black	M-N318R19	12 pcs	\$22.00
		4/0	18" (45cm)	Black	M-N418R19	12 pcs	\$22.00
		5/0	18" (45cm)	Black	M-N518R19	12 pcs	\$22.00
		4/0	10" (25cm)	Black	M-N410R19-B	12 pcs	\$21.00
		5/0	10" (25cm)	Black	M-N510R19-B	12 pcs	\$21.00
		3/0	18" (45cm)	Black	M-N318XR19	36 pcs	\$63.00
		4/0	18" (45cm)	Black	M-N418XR19	36 pcs	\$63.00
		5/0	18" (45cm)	Black	M-N518XR19	36 pcs	\$63.00
<b>MEDIUM</b> FS-3 / C-12 / C22 16mm, 3/8 Circle		4/0	10" (25cm)	Black	M-N410R13-B	12 pcs	\$21.00
		5/0	10" (25cm)	Black	S-N510R13-B	12 pcs	\$21.00
		6/0	10" (25cm)	Black	S-N610R13-B	12 pcs	\$21.00
		4/0	18" (45cm)	Black	M-N418R16	12 pcs	\$22.00
<b>SMALL</b> P-3 / P-13 / C3 13mm, 3/8 Circle		5/0	18" (45cm)	Black	M-N518R16	12 pcs	\$22.00
		6/0	18" (45cm)	Black	M-N618R13	12 pcs	\$22.00
		4/0	18" (45cm)	Black	S-N418R13	12 pcs	\$22.00
		5/0	18" (45cm)	Black	S-N518R13	12 pcs	\$22.00
		6/0	18" (45cm)	Black	S-N618R13	12 pcs	\$22.00
		4/0	18" (45cm)	Black	S-N418XR13	36 pcs	\$63.00
<b>X-SMALL</b> P-1 / P-10 / C1 11mm, 3/8 Circle		5/0	18" (45cm)	Black	S-N518XR13	36 pcs	\$63.00
		6/0	18" (45cm)	Black	S-N618XR13	36 pcs	\$63.00
		6/0	18" (45cm)	Black	XS-N618R11	12 pcs	\$41.50



**POLYESTER** Comparable to Ethibond® \* and Surgidac® \*\*

Needle Size	Needle Shape	Thread Size (USP)	Thread Length	Thread Color	Item #	Qty.	Price/box
<b>X-LARGE</b> FSL / C-16 30mm, 3/8 Circle		0	30" (75cm)	Green	XL-E030R30	12 pcs	\$38.00
<b>LARGE</b> FS-1 / C-14 / C7 24mm, 3/8 Circle		2/0	30" (75cm)	Green	L-E230R24	12 pcs	\$38.00
<b>MEDIUM</b> FS-2 / C-13 / C6 19mm, 3/8 Circle		3/0	18" (45cm)	Green	M-E318R19	12 pcs	\$32.00






**PREMIUM POLYPROPYLENE** Comparable to Prolene® \* and Surgipro II® \*\*

<b>LARGE</b> PS-1 / P-14 / PC33 24mm, 3/8 Circle		3/0 4/0	27" (69cm) 27" (69cm)	Blue Blue	PLP-327R24 PLP-427R24	12 pcs 12 pcs	\$55.50 \$55.50
<b>MEDIUM</b> PS-2 / P-12 / PC31 19mm, 3/8 Circle		3/0 4/0 5/0	18" (45cm) 18" (45cm) 18" (45cm)	Blue Blue Blue	PMP-318R19 PMP-418R19 PMP-518R19	12 pcs 12 pcs 12 pcs	\$49.00 \$49.00 \$49.00
<b>MEDIUM</b> PS-3 / P-11 / PC34 16mm, 3/8 Circle		4/0 5/0	18" (45cm) 18" (45cm)	Blue Blue	PMP-418R16 PMP-518R16	12 pcs 12 pcs	\$49.00 \$49.00
<b>SMALL</b> P-3 / P-13 / C3 13mm, 3/8 Circle		4/0 5/0 6/0	18" (45cm) 18" (45cm) 18" (45cm)	Blue Blue Blue	PSP-418R13 PSP-518R13 PSP-618R13	12 pcs 12 pcs 12 pcs	\$49.00 \$49.00 \$49.00






**POLYPROPYLENE** Comparable to Prolene® \* and Surgipro II® \*\*





<b>X-LARGE</b> CT-1 / GS-21 / T8 36mm, 1/2 Circle		0	30" (75cm)	Blue	XL-P030T36	12 pcs	\$38.00
<b>LARGE</b> FS-1 / C-14 / C7 24mm, 3/8 Circle		2/0	30" (75cm)	Blue	L-P230R24	12 pcs	\$38.00


**POLYPROPYLENE** Comparable to Prolene® \* and Surgipro II® \*\*




<b>MEDIUM</b> FS-2 / C-13 / C6 19mm, 3/8 Circle		3/0 4/0	30" (75cm) 30" (75cm)	Blue Blue	M-P330R19 M-P430R19	12 pcs 12 pcs	\$38.00 \$38.00
<b>MEDIUM</b> FS-3 / C-12 / C22 16mm, 3/8 Circle		3/0 4/0 5/0 4/0	18" (45cm) 18" (45cm) 18" (45cm) 18" (45cm)	Blue Blue Blue Blue	M-P318R19 M-P418R19 M-P518R19 M-P418XR19	12 pcs 12 pcs 12 pcs 36 pcs	\$31.00 \$31.00 \$31.00 \$90.00
<b>MEDIUM</b> FS-3 / C-12 / C22 16mm, 3/8 Circle		4/0 5/0	18" (45cm) 18" (45cm)	Blue Blue	M-P418R16 M-P518R16	12 pcs 12 pcs	\$31.00 \$31.00
<b>SMALL</b> P-3 / P-13 / C3 13mm, 3/8 Circle		4/0 5/0 6/0 4/0 5/0 6/0	18" (45cm) 18" (45cm) 18" (45cm) 18" (45cm) 18" (45cm) 18" (45cm)	Blue Blue Blue Blue Blue Blue	S-P418R13 S-P518R13 S-P618R13 S-P418XR13 S-P518XR13 S-P618XR13	12 pcs 12 pcs 12 pcs 36 pcs 36 pcs 36 pcs	\$31.00 \$31.00 \$31.00 \$90.00 \$90.00 \$90.00
<b>X-SMALL</b> P-1 / P-10 / C1 11mm, 3/8 Circle		6/0 7/0	18" (45cm) 18" (45cm)	Blue Blue	XS-P618R11 XS-P718R11	12 pcs 12 pcs	\$41.50 \$41.50



**PTFE** Comparable to Cytoplast® \*\*\* PTFE sutures



<b>XX-LARGE</b> KS / SC-2 60mm		2/0	30" (75cm)	White	XXL-F230K60	12 pcs	\$78.00
<b>X-LARGE</b> CT-1 / GS-21 / T8 36mm, 1/2 Circle		2/0	30" (75cm)	White	XL-F230T36	12 pcs	\$78.00
<b>MEDIUM</b> FS-2 / C-13 / C6 19mm, 3/8 Circle		2/0 3/0 4/0	18" (45cm) 18" (45cm) 18" (45cm)	White White White	M-F218R19 M-F318R19 M-F418R19	12 pcs 12 pcs 12 pcs	\$78.00 \$78.00 \$78.00
<b>MEDIUM</b> FS-3 / C-12 / C22 16mm, 3/8 Circle		3/0 4/0 5/0	18" (45cm) 18" (45cm) 18" (45cm)	White White White	M-F318R16 M-F418R16 M-F518R16	12 pcs 12 pcs 12 pcs	\$78.00 \$78.00 \$78.00
<b>SMALL</b> P-3 / P-13 / C3 13mm, 3/8 Circle		4/0 5/0 6/0	18" (45cm) 18" (45cm) 18" (45cm)	White White White	S-F418R13 S-F518R13 S-F618R13	12 pcs 12 pcs 12 pcs	\$78.00 \$78.00 \$78.00









Needle Size	Needle Shape	Thread Size (USP)	Thread Length	Thread Color	Item #	Qty.	Price/box
<b>SILK</b> <span style="float: right;">Comparable to Perma-hand® * and Sofsilk® **</span>							
<b>XX-LARGE</b> KS / SC-2 60mm		2/0	30" (75cm)	Black	XXL-S230K60	12 pcs	\$25.00
		3/0	30" (75cm)	Black	M-S330R19	12 pcs	\$25.00
		4/0	30" (75cm)	Black	M-S430R19	12 pcs	\$25.00
<b>MEDIUM</b> FS-2 / C-13 / C6 19mm, 3/8 Circle		3/0	18" (45cm)	Black	M-S318R19	12 pcs	\$21.00
		4/0	18" (45cm)	Black	M-S418R19	12 pcs	\$21.00
		5/0	18" (45cm)	Black	M-S518R19	12 pcs	\$21.00
<b>MEDIUM</b> FS-3 / C-12 / C22 16mm, 3/8 Circle		4/0	18" (45cm)	Black	M-S418R16	12 pcs	\$21.00
		5/0	18" (45cm)	Black	M-S518R16	12 pcs	\$21.00
<b>SMALL</b> P-3 / P-13 / C3 13mm, 3/8 Circle		4/0	18" (45cm)	Black	S-S418R13	12 pcs	\$21.00
		5/0	18" (45cm)	Black	S-S518R13	12 pcs	\$21.00
		6/0	18" (45cm)	Black	S-S618R13	12 pcs	\$21.00








<b>SURGICAL STAINLESS STEEL</b>							
<b>XX-LARGE</b> CTX / GS-25 48mm, 1/2 Circle		4	18" (45cm)	S.Steel	W4D-18T48	12 pcs	\$99.50
		6	18" (45cm)	S.Steel	W6Q-18T48	12 pcs	\$169.50

<b>PREMIUM PDO</b> <span style="float: right;">Comparable to PDS II® * and Maxon® **</span>								
<b>LARGE</b> PS-1 / P-14 / PC33 24mm, 3/8 Circle		2/0	27" (69cm)	Violet	PLD-227R24	12 pcs	\$59.50	
		3/0	27" (69cm)	Violet	PLD-327R24	12 pcs	\$59.50	
		4/0	27" (69cm)	Violet	PLD-427R24	12 pcs	\$59.50	
<b>MEDIUM</b> PS-2 / P-12 / PC31 19mm, 3/8 Circle		2/0	27" (69cm)	Violet	PMD-227R19	12 pcs	\$59.50	
		3/0	27" (69cm)	Violet	PMD-327R19	12 pcs	\$59.50	
		4/0	27" (69cm)	Violet	PMD-427R19	12 pcs	\$59.50	
		2/0	18" (45cm)	Violet	PMD-218R19	12 pcs	\$50.50	
		3/0	18" (45cm)	Violet	PMD-318R19	12 pcs	\$50.50	
<b>MEDIUM</b> PS-3 / P-11 / PC34 16mm, 3/8 Circle		4/0	18" (45cm)	Violet	PMD-418R16	12 pcs	\$50.50	
		<b>SMALL</b> P-3 / P-13 / C3 13mm, 3/8 Circle	4/0	18" (45cm)	Violet	PSD-418R13	12 pcs	\$50.50
			5/0	18" (45cm)	Violet	PSD-518R13	12 pcs	\$50.50





Needle Size	Needle Shape	Thread Size (USP)	Thread Length	Thread Color	Item #	Qty.	Price/box	
<b>PDO</b> <span style="float: right;">Comparable to PDS II® * and Maxon® **</span>								
<b>XX-LARGE</b> CP / GS-12 / C14 40mm, 1/2 Circle		2	30" (75cm)	Violet	XXL-2D30R40	12 pcs	\$51.00	
		1	30" (75cm)	Violet	XXL-1D30R40	12 pcs	\$51.00	
<b>X-LARGE</b> CT-1 / GS-21 / T8 36mm, 1/2 Circle		0	30" (75cm)	Violet	XL-D030T36	12 pcs	\$51.00	
		2/0	30" (75cm)	Violet	XL-D230T36	12 pcs	\$51.00	
<b>X-LARGE</b> FSL / C-16 30mm, 3/8 Circle		0	30" (75cm)	Violet	XL-D030R30	12 pcs	\$51.00	
		<b>LARGE</b> SH (V-20) 26mm, 1/2 Circle	2/0	30" (75cm)	Violet	L-D230T26	12 pcs	\$47.00
3/0	30" (75cm)		Violet	L-D330T26	12 pcs	\$47.00		
<b>LARGE</b> FS-1 / C-14 / C7 24mm, 3/8 Circle		2/0	30" (75cm)	Violet	L-D230R24	12 pcs	\$47.00	
		3/0	30" (75cm)	Violet	L-D330R24	12 pcs	\$47.00	
		2/0	30" (75cm)	Violet	L-D230XR24	36 pcs	\$133.00	
		3/0	30" (75cm)	Violet	L-D330XR24	36 pcs	\$133.00	
		<b>MEDIUM</b> FS-2 / C-13 / C6 19mm, 3/8 Circle	2/0	30" (75cm)	Violet	M-D230R19	12 pcs	\$47.00
			3/0	30" (75cm)	Violet	M-D330R19	12 pcs	\$47.00
4/0	30" (75cm)		Violet	M-D430R19	12 pcs	\$47.00		
3/0	18" (45cm)		Violet	M-D318R19	12 pcs	\$42.00		
4/0	18" (45cm)		Violet	M-D418R19	12 pcs	\$42.00		
<b>MEDIUM</b> RB-1 / CV-23 / T28 17mm, 1/2 Circle		5/0	18" (45cm)	Violet	M-D518R19	12 pcs	\$42.00	
		3/0	30" (75cm)	Violet	M-D330T17	12 pcs	\$47.00	
<b>MEDIUM</b> FS-3 / C-12 / C22 16mm, 3/8 Circle		4/0	18" (45cm)	Violet	M-D418R16	12 pcs	\$42.00	
		5/0	18" (45cm)	Violet	M-D518R16	12 pcs	\$42.00	
<b>SMALL</b> P-3 / P-13 / C3 13mm, 3/8 Circle		4/0	18" (45cm)	Violet	S-D418R13	12 pcs	\$42.00	
		5/0	18" (45cm)	Violet	S-D518R13	12 pcs	\$42.00	
		6/0	18" (45cm)	Violet	S-D618R13	12 pcs	\$42.00	






<b>Knotless Barbed PDO</b> <span style="float: right;">Comparable to Stratafix PDO® * and V-Loc® **</span>							
<b>X-LARGE</b> CT-1 / GS-21 / T8 36mm, 1/2 Circle		2/0	18" (45cm)	Violet	B-D218T36	6 pcs	\$139.50
		<b>LARGE</b> SH (V-20) 26mm, 1/2 Circle	2/0	18" (45cm)	Violet	B-D218T26	6 pcs
<b>MEDIUM</b> PS-2 / P-12 / PC31 19mm, 3/8 Circle			3/0	12" (30cm)	Violet	B-PD312R19	6 pcs

Needle Size	Needle Shape	Thread Size (USP)	Thread Length	Thread Color	Item #	Qty.	Price/box
<b>PGA MICROSUTURES</b> <span style="float: right;">Comparable to Vicryl® *and Polysorb® **</span>							
<b>MICRO</b> TG-140-8 (SE-140-6) 6.4mm, 3/8 Circle		6/0	12" (30cm)	Violet	XXS-G612S6D	12 pcs	\$169.50
		8/0	12" (30cm)	Violet	XXS-G812S6D	12 pcs	\$169.50
<b>PREMIUM PGA</b> <span style="float: right;">Comparable to Vicryl® *and Polysorb® **</span>							
<b>XX-LARGE</b> CT / GS-24 / T10 40mm, 1/2 Circle		2/0	27" (69cm)	Violet	PXXL-G227T40	12 pcs	\$59.50
<b>X-LARGE</b> CT-1 / GS-21 / T8 36mm, 1/2 Circle		2/0	27" (69cm)	Undyed	PXLGU-227T36	12 pcs	\$59.50
<b>LARGE</b> SH (V-20) 26mm, 1/2 Circle		3/0	27" (69cm)	Undyed	PLGU-327T26	12 pcs	\$59.50
<b>LARGE</b> PS-1 / P-14 / PC33 24mm, 3/8 Circle		3/0	27" (69cm)	Violet	PLG-327R24	12 pcs	\$59.50
		4/0	27" (69cm)	Violet	PLG-427R24	12 pcs	\$59.50
		2/0	27" (69cm)	Undyed	PLGU-227R24	12 pcs	\$59.50
		3/0	27" (69cm)	Undyed	PLGU-327R24	12 pcs	\$59.50
<b>MEDIUM</b> PS-2 / P-12 / PC31 19mm, 3/8 Circle		4/0	27" (69cm)	Violet	PMG-427R19	12 pcs	\$59.50
		4/0	27" (69cm)	Undyed	PMGU-427R19	12 pcs	\$59.50
		3/0	18" (45cm)	Violet	PMG-318R19	12 pcs	\$50.50
		4/0	18" (45cm)	Violet	PMG-418R19	12 pcs	\$50.50
		5/0	18" (45cm)	Violet	PMG-518R19	12 pcs	\$50.50
		3/0	18" (45cm)	Undyed	PMGU-318R19	12 pcs	\$50.50
		4/0	18" (45cm)	Undyed	PMGU-418R19	12 pcs	\$50.50
		5/0	18" (45cm)	Undyed	PMGU-518R19	12 pcs	\$50.50
<b>MEDIUM</b> PS-3 / P-11 / PC34 16mm, 3/8 Circle		4/0	18" (45cm)	Undyed	PMGU-418R16	12 pcs	\$50.50
		5/0	18" (45cm)	Undyed	PMGU-518R16	12 pcs	\$50.50
		4/0	18" (45cm)	Violet	PMG-418R16	12 pcs	\$50.50
		5/0	18" (45cm)	Violet	PMG-518R16	12 pcs	\$50.50
<b>SMALL</b> P-3 / P-13 / C3 13mm, 3/8 Circle		4/0	18" (45cm)	Undyed	PSGU-418R13	12 pcs	\$50.50
		5/0	18" (45cm)	Undyed	PSGU-518R13	12 pcs	\$50.50
		6/0	18" (45cm)	Undyed	PSGU-618R13	12 pcs	\$50.50
		4/0	18" (45cm)	Violet	PSG-418R13	12 pcs	\$50.50
		5/0	18" (45cm)	Violet	PSG-518R13	12 pcs	\$50.50
		6/0	18" (45cm)	Violet	PSG-618R13	12 pcs	\$50.50

Needle Size	Needle Shape	Thread Size (USP)	Thread Length	Thread Color	Item #	Qty.	Price/box
<b>PGA</b> <span style="float: right;">Comparable to Vicryl® *and Polysorb® **</span>							
<b>XX-LARGE</b> CT / GS-24 / T10 40mm, 1/2 Circle		1	30" (75cm)	Violet	XXL-1G30T40	12 pcs	\$51.00
		0	30" (75cm)	Violet	XXL-G030T40	12 pcs	\$51.00
		2/0	30" (75cm)	Violet	XXL-G230T40	12 pcs	\$51.00
		1	30" (75cm)	Undyed	XXL-1G30T40-U	12 pcs	\$51.00
		0	30" (75cm)	Undyed	XXL-G030T40-U	12 pcs	\$51.00
2/0	30" (75cm)	Undyed	XXL-G230T40-U	12 pcs	\$51.00		
<b>X-LARGE</b> CT-1 / GS-21 / T8 36mm, 1/2 Circle		0	30" (75cm)	Violet	XL-G030T36	12 pcs	\$51.00
		2/0	30" (75cm)	Undyed	XL-G230T36-U	12 pcs	\$51.00
<b>X-LARGE</b> FSL / C-16 30mm, 3/8 Circle		0	30" (75cm)	Violet	XL-G030R30	12 pcs	\$51.00
<b>LARGE</b> SH (V-20) 26mm, 1/2 Circle		2/0	30" (75cm)	Violet	L-G230T26	12 pcs	\$47.00
		2/0	30" (75cm)	Undyed	L-G230T26-U	12 pcs	\$47.00
		3/0	30" (75cm)	Violet	L-G330T26	12 pcs	\$47.00
<b>LARGE</b> FS-1 / C-14 / C7 24mm, 3/8 Circle		2/0	30" (75cm)	Violet	L-G230R24	12 pcs	\$47.00
		3/0	30" (75cm)	Violet	L-G330R24	12 pcs	\$47.00
		2/0	30" (75cm)	Undyed	L-G230R24-U	12 pcs	\$47.00
		3/0	30" (75cm)	Undyed	L-G330R24-U	12 pcs	\$47.00
		3/0	30" (75cm)	Violet	M-G330R19	12 pcs	\$47.00
		4/0	30" (75cm)	Violet	M-G430R19	12 pcs	\$47.00
<b>MEDIUM</b> FS-2 / C-13 / C6 19mm, 3/8 Circle		3/0	30" (75cm)	Undyed	M-G330R19-U	12 pcs	\$47.00
		4/0	30" (75cm)	Undyed	M-G430R19-U	12 pcs	\$47.00
		3/0	18" (45cm)	Violet	M-G318R19	12 pcs	\$42.00
		4/0	18" (45cm)	Violet	M-G418R19	12 pcs	\$42.00
		5/0	18" (45cm)	Violet	M-G518R19	12 pcs	\$42.00
		3/0	18" (45cm)	Undyed	M-G318R19-U	12 pcs	\$42.00
		4/0	18" (45cm)	Undyed	M-G418R19-U	12 pcs	\$42.00
		5/0	18" (45cm)	Undyed	M-G518R19-U	12 pcs	\$42.00
<b>MEDIUM</b> RB-1 / CV-23 / T28 17mm, 1/2 Circle		3/0	30" (75cm)	Undyed	M-G330T17-U	12 pcs	\$47.00
		4/0	30" (75cm)	Violet	M-G430T17	12 pcs	\$47.00
		4/0	30" (75cm)	Undyed	M-G430T17-U	12 pcs	\$47.00

Needle Size	Needle Shape	Thread Size (USP)	Thread Length	Thread Color	Item #	Qty.	Price/box
<b>PGA</b> <span style="float: right;">Comparable to Vicryl® *and Polysorb® **</span>							
<b>MEDIUM</b> FS-3 / C-12 / C22 16mm, 3/8 Circle		4/0	18" (45cm)	Violet	M-G418R16	12 pcs	\$42.00
		5/0	18" (45cm)	Violet	M-G518R16	12 pcs	\$42.00
		4/0	18" (45cm)	Undyed	M-G418R16-U	12 pcs	\$42.00
		5/0	18" (45cm)	Undyed	M-G518R16-U	12 pcs	\$42.00
		<b>MEDIUM</b> PC-3 / PC-11 16mm, 3/8 Circle					
		4/0	18" (45cm)	Undyed	M-G418A16-U	12 pcs	\$42.00
		4/0	18" (45cm)	Violet	S-G418R13	12 pcs	\$42.00
		5/0	18" (45cm)	Violet	S-G518R13	12 pcs	\$42.00
		6/0	18" (45cm)	Violet	S-G618R13	12 pcs	\$42.00
		<b>SMALL</b> P-3 / P-13 / C3 13mm, 3/8 Circle					
		4/0	18" (45cm)	Undyed	S-G418R13-U	12 pcs	\$42.00
		5/0	18" (45cm)	Undyed	S-G518R13-U	12 pcs	\$42.00
		6/0	18" (45cm)	Undyed	S-G618R13-U	12 pcs	\$42.00
		4/0	18" (45cm)	Violet	S-G418XR13	36 pcs	\$119.00
		4/0	18" (45cm)	Undyed	S-G418XR13-U	36 pcs	\$119.00
5/0	18" (45cm)	Undyed	S-G518XR13-U	36 pcs	\$119.00		
<b>X-SMALL</b> P-1 / P-10 / C1 11mm, 3/8 Circle							
		5/0	18" (45cm)	Undyed	XS-G518R11-U	12 pcs	\$51.50
		6/0	18" (45cm)	Undyed	XS-G618R11-U	12 pcs	\$51.50
		7/0	18" (45cm)	Violet	XS-G718R11	12 pcs	\$51.50
<b>PREMIUM RAPID PGA</b> <span style="float: right;">Comparable to Vicryl Rapide* and Velosorb Fast**</span>							
<b>MEDIUM</b> PS-2 / P-12 / PC31 19mm, 3/8 Circle							
		4/0	18" (45cm)	Undyed	PMK-418R19	12 pcs	\$52.00
		<b>SMALL</b> P-3 / P-13 / C3 13mm, 3/8 Circle					
		4/0	18" (45cm)	Undyed	PSK-418R13	12 pcs	\$52.00
		5/0	18" (45cm)	Undyed	PSK-518R13	12 pcs	\$52.00
<b>RAPID PGA</b> <span style="float: right;">Comparable to Vicryl Rapide* and Velosorb Fast**</span>							
<b>MEDIUM</b> FS-2 / C-13 / C6 19mm, 3/8 Circle							
		4/0	30" (75cm)	Undyed	M-K430R19	12 pcs	\$51.00
		3/0	18" (45cm)	Undyed	M-K318R19	12 pcs	\$42.00
		4/0	18" (45cm)	Undyed	M-K418R19	12 pcs	\$42.00
		5/0	18" (45cm)	Undyed	M-K518R19	12 pcs	\$42.00
<b>SMALL</b> P-3 / P-13 / C3 13mm, 3/8 Circle							
		4/0	18" (45cm)	Undyed	S-K418R13	12 pcs	\$42.00
		5/0	18" (45cm)	Undyed	S-K518R13	12 pcs	\$42.00
<b>X-SMALL</b> P-1 / P-10 / C1 11mm, 3/8 Circle							
		5/0	18" (45cm)	Undyed	XS-K518R11	12 pcs	\$51.50











Needle Size	Needle Shape	Thread Size (USP)	Thread Length	Thread Color	Item #	Qty.	Price/box
<b>PREMIUM PGCL</b> <span style="float: right;">Comparable to Monocryl® * and Biosyn®**</span>							
<b>LARGE</b> PS-1 / P-14 / PC33 24mm, 3/8 Circle							
		3/0	27" (69cm)	Violet	PLQC-327R24	12 pcs	\$59.50
		4/0	27" (69cm)	Violet	PLQC-427R24	12 pcs	\$59.50
		2/0	27" (69cm)	Undyed	PLQ-227R24	12 pcs	\$59.50
		3/0	27" (69cm)	Undyed	PLQ-327R24	12 pcs	\$59.50
		4/0	27" (69cm)	Undyed	PLQ-427R24	12 pcs	\$59.50
<b>MEDIUM</b> PS-2 / P-12 / PC31 19mm, 3/8 Circle							
		3/0	27" (69cm)	Undyed	PMQ-327R19	12 pcs	\$59.50
		4/0	27" (69cm)	Undyed	PMQ-427R19	12 pcs	\$59.50
		3/0	18" (45cm)	Violet	PMQC-318R19	12 pcs	\$50.50
		4/0	18" (45cm)	Violet	PMQC-418R19	12 pcs	\$50.50
		5/0	18" (45cm)	Violet	PMQC-518R19	12 pcs	\$50.50
<b>SMALL</b> P-3 / P-13 / C3 13mm, 3/8 Circle							
		3/0	18" (45cm)	Undyed	PMQ-318R19	12 pcs	\$50.50
		4/0	18" (45cm)	Undyed	PMQ-418R19	12 pcs	\$50.50
		5/0	18" (45cm)	Undyed	PMQ-518R19	12 pcs	\$50.50
		4/0	18" (45cm)	Violet	PSQC-418R13	12 pcs	\$50.50
		5/0	18" (45cm)	Violet	PSQC-518R13	12 pcs	\$50.50
<b>SMALL</b> TF / CV-22 / T22 13mm, 1/2 Circle							
		4/0	18" (45cm)	Undyed	PSQ-418R13	12 pcs	\$50.50
		5/0	18" (45cm)	Undyed	PSQ-518R13	12 pcs	\$50.50
		6/0	18" (45cm)	Undyed	PSQ-618R13	12 pcs	\$50.50
		5/0	18" (45cm)	Undyed	PSQ-518T13	12 pcs	\$50.50

Needle Size	Needle Shape	Thread Size (USP)	Thread Length	Thread Color	Item #	Qty.	Price/box
<b>PGCL</b> <span style="float: right;">Comparable to Monocryl® * and Biosyn®**</span>							
<b>X-LARGE</b> CT-1 / GS-21 / T8 36mm, 1/2 Circle		2/0	30" (75cm)	Undyed	XL-Q230T36	12 pcs	\$51.00
<b>LARGE</b> FS-1 / C-14 / C7 24mm, 3/8 Circle		2/0	30" (75cm)	Undyed	L-Q230R24	12 pcs	\$51.00
		3/0	30" (75cm)	Undyed	L-Q330R24	12 pcs	\$51.00
		4/0	30" (75cm)	Undyed	L-Q430R24	12 pcs	\$51.00
		3/0	30" (75cm)	Undyed	M-Q330R19	12 pcs	\$51.00
		4/0	30" (75cm)	Undyed	M-Q430R19	12 pcs	\$51.00
		3/0	18" (45cm)	Undyed	M-Q318R19	12 pcs	\$42.00
		4/0	18" (45cm)	Undyed	M-Q418R19	12 pcs	\$42.00
<b>MEDIUM</b> FS-2 / C-13 / C6 19mm, 3/8 Circle		5/0	18" (45cm)	Undyed	M-Q518R19	12 pcs	\$42.00
		3/0	18" (45cm)	Violet	M-Q318R19-C	12 pcs	\$42.00
		4/0	18" (45cm)	Violet	M-Q418R19-C	12 pcs	\$42.00
		5/0	18" (45cm)	Violet	M-Q518R19-C	12 pcs	\$42.00
		3/0	30" (75cm)	Undyed	M-Q330XR19	36 pcs	\$145.00
		4/0	30" (75cm)	Undyed	M-Q430XR19	36 pcs	\$145.00
		3/0	18" (45cm)	Undyed	M-Q318XR19	36 pcs	\$119.00
<b>SMALL</b> P-3 / P-13 / C3 13mm, 3/8 Circle		4/0	18" (45cm)	Undyed	M-Q418XR19	36 pcs	\$119.00
		4/0	18" (45cm)	Violet	M-Q418XR19-C	36 pcs	\$119.00
		4/0	18" (45cm)	Undyed	S-Q418R13	12 pcs	\$42.00
		5/0	18" (45cm)	Undyed	S-Q518R13	12 pcs	\$42.00
		6/0	18" (45cm)	Undyed	S-Q618R13	12 pcs	\$42.00
		4/0	18" (45cm)	Violet	S-Q418R13-C	12 pcs	\$42.00
		5/0	18" (45cm)	Violet	S-Q518R13-C	12 pcs	\$42.00
		6/0	18" (45cm)	Violet	S-Q618R13-C	12 pcs	\$42.00
		4/0	18" (45cm)	Undyed	S-Q418XR13	36 pcs	\$119.00
		5/0	18" (45cm)	Undyed	S-Q518XR13	36 pcs	\$119.00
<b>X-SMALL</b> P-1 / P-10 / C1 11mm, 3/8 Circle		5/0	18" (45cm)	Violet	XS-Q518R11-C	12 pcs	\$51.50











\* Registered trademarks of Ethicon    \*\* Registered trademarks of Covidien    \*\*\* Registered trademark of Osteogenics Biomedical

Needle Size	Needle Shape	Thread Size (USP)	Thread Length	Thread Color	Item #	Qty.	Price/box
<b>PREMIUM CHROMIC GUT</b>							
<b>LARGE</b> PS-1 / P-14 / PC33 24mm, 3/8 Circle		2/0	27" (69cm)	Brown	PLC-227R24	12 pcs	\$49.50
		3/0	27" (69cm)	Brown	PLC-327R24	12 pcs	\$49.50
<b>MEDIUM</b> PS-2 / P-12 / PC31 19mm, 3/8 Circle		3/0	27" (69cm)	Brown	PMC-327R19	12 pcs	\$49.50
		4/0	27" (69cm)	Brown	PMC-427R19	12 pcs	\$49.50
		3/0	18" (45cm)	Brown	PMC-318R19	12 pcs	\$48.00
		4/0	18" (45cm)	Brown	PMC-418R19	12 pcs	\$48.00
<b>MEDIUM</b> PS-3 / P-11 / PC34 16mm, 3/8 Circle		4/0	18" (45cm)	Brown	PMC-418R16	12 pcs	\$48.00
		5/0	18" (45cm)	Brown	PMC-518R16	12 pcs	\$48.00
<b>SMALL</b> P-3 / P-13 / C3 13mm, 3/8 Circle		4/0	18" (45cm)	Brown	PSC-418R13	12 pcs	\$48.00
		5/0	18" (45cm)	Brown	PSC-518R13	12 pcs	\$48.00
		6/0	18" (45cm)	Brown	PSC-618R13	12 pcs	\$48.00
<b>CHROMIC GUT</b>							
<b>LARGE</b> FS-1 / C-14 / C7 24mm, 3/8 Circle		2/0	30" (75cm)	Brown	L-C230R24	12 pcs	\$39.50
<b>MEDIUM</b> FS-2 / C-13 / C6 19mm, 3/8 Circle		3/0	30" (75cm)	Brown	M-C330R19	12 pcs	\$39.50
		4/0	30" (75cm)	Brown	M-C430R19	12 pcs	\$39.50
		3/0	18" (45cm)	Brown	M-C318R19	12 pcs	\$38.00
		4/0	18" (45cm)	Brown	M-C418R19	12 pcs	\$38.00
		5/0	18" (45cm)	Brown	M-C518R19	12 pcs	\$38.00
<b>MEDIUM</b> FS-3 / C-12 / C22 16mm, 3/8 Circle		4/0	18" (45cm)	Brown	M-C418R16	12 pcs	\$38.00
		5/0	18" (45cm)	Brown	M-C518R16	12 pcs	\$38.00
<b>SMALL</b> P-3 / P-13 / C3 13mm, 3/8 Circle		4/0	18" (45cm)	Brown	S-C418R13	12 pcs	\$38.00
		5/0	18" (45cm)	Brown	S-C518R13	12 pcs	\$38.00
		6/0	18" (45cm)	Brown	S-C618R13	12 pcs	\$38.00



Needle Size	Needle Shape	Thread Size (USP)	Thread Length	Thread Color	Item #	Qty.	Price/box
<b>PREMIUM FAST GUT</b>							
SMALL P-3 / P-13 / C3 13mm, 3/8 Circle		5/0	18" (45cm)	Beige	PSRJ-518R13	12 pcs	\$49.50
		6/0	18" (45cm)	Beige	PSRJ-618R13	12 pcs	\$49.50
SMALL PC-1 13mm, 3/8 Circle		5/0	18" (45cm)	Beige	PSRJ-518A13	12 pcs	\$49.50
<b>PREMIUM PLAIN GUT</b>							
LARGE PS-1 / P-14 / PC33 24mm, 3/8 Circle		3/0	27" (69cm)	Beige	PLJ-327R24	12 pcs	\$49.50
MEDIUM PS-2 / P-12 / PC31 19mm, 3/8 Circle		3/0	27" (69cm)	Beige	PMJ-327R19	12 pcs	\$49.50
		4/0	27" (69cm)	Beige	PMJ-427R19	12 pcs	\$49.50
		3/0	18" (45cm)	Beige	PMJ-318R19	12 pcs	\$48.00
		4/0	18" (45cm)	Beige	PMJ-418R19	12 pcs	\$48.00
MEDIUM PS-3 / P-11 / PC34 16mm, 3/8 Circle		4/0	18" (45cm)	Beige	PMJ-418R16	12 pcs	\$48.00
		5/0	18" (45cm)	Beige	PMJ-518R16	12 pcs	\$48.00
SMALL P-3 / P-13 / C3 13mm, 3/8 Circle		4/0	18" (45cm)	Beige	PSJ-418R13	12 pcs	\$48.00
		5/0	18" (45cm)	Beige	PSJ-518R13	12 pcs	\$48.00
		6/0	18" (45cm)	Beige	PSJ-618R13	12 pcs	\$48.00
<b>PLAIN GUT</b>							
LARGE FS-1 / C-14 / C7 24mm, 3/8 Circle		2/0	30" (75cm)	Beige	L-J230R24	12 pcs	\$39.50
MEDIUM FS-2 / C-13 / C6 19mm, 3/8 Circle		3/0	30" (75cm)	Beige	M-J330R19	12 pcs	\$39.50
		4/0	30" (75cm)	Beige	M-J430R19	12 pcs	\$39.50
		3/0	18" (45cm)	Beige	M-J318R19	12 pcs	\$38.00
		4/0	18" (45cm)	Beige	M-J418R19	12 pcs	\$38.00
		5/0	18" (45cm)	Beige	M-J518R19	12 pcs	\$38.00
MEDIUM FS-3 / C-12 / C22 16mm, 3/8 Circle		4/0	18" (45cm)	Beige	M-J418R16	12 pcs	\$38.00
		5/0	18" (45cm)	Beige	M-J518R16	12 pcs	\$38.00
SMALL P-3 / P-13 / C3 13mm, 3/8 Circle		4/0	18" (45cm)	Beige	S-J418R13	12 pcs	\$38.00
		5/0	18" (45cm)	Beige	S-J518R13	12 pcs	\$38.00
		6/0	18" (45cm)	Beige	S-J618R13	12 pcs	\$38.00

## UNIFY® VET SUTURES For Veterinary Use Only

Needle Size	Needle Shape	Thread Size (USP)	Thread Length	Thread Color	Item #	Qty.	Price/box
<b>VET NYLON</b> <span style="float: right;">Comparable to Ethilon® * and Monosof® **</span>							
X-LARGE CP-1 / GS-11 36mm, 1/2 Circle		0	30" (75cm)	Black	VETN-030R36	12 pcs	\$20.00
		0	30" (75cm)	Black	VETN-030R30	12 pcs	\$20.00
X-LARGE FSL / C-16 30mm, 3/8 Circle		0	30" (75cm)	Black	VETN-030R30	12 pcs	\$20.00
		2/0	30" (75cm)	Black	VETN-230R24	12 pcs	\$20.00
LARGE FS-1 / C-14 / C7 24mm, 3/8 Circle		3/0	30" (75cm)	Black	VETN-330R24	12 pcs	\$20.00
		4/0	30" (75cm)	Black	VETN-430R19	12 pcs	\$20.00
MEDIUM FS-2 / C-13 / C6 19mm, 3/8 Circle		4/0	30" (75cm)	Black	VETN-430R19	12 pcs	\$20.00
<b>VET PGA</b> <span style="float: right;">Comparable to Vicryl® * and Polysorb® **</span>							
X-LARGE CP-1 / GS-11 36mm, 1/2 Circle		0	30" (75cm)	Violet	VETG-030R36	12 pcs	\$32.00
		0	30" (75cm)	Violet	VETG-030R30	12 pcs	\$32.00
X-LARGE FSL / C-16 30mm, 3/8 Circle		0	30" (75cm)	Violet	VETG-030R30	12 pcs	\$32.00
		2/0	30" (75cm)	Violet	VETG-230R24	12 pcs	\$32.00
LARGE FS-1 / C-14 / C7 24mm, 3/8 Circle		3/0	30" (75cm)	Violet	VETG-330R24	12 pcs	\$32.00
		4/0	30" (75cm)	Violet	VETG-430R19	12 pcs	\$32.00
MEDIUM FS-2 / C-13 / C6 19mm, 3/8 Circle		4/0	30" (75cm)	Violet	VETG-430R19	12 pcs	\$32.00
<b>VET PDO</b> <span style="float: right;">Comparable to PDS II® * and Maxon® **</span>							
X-LARGE CP-1 / GS-11 36mm, 1/2 Circle		0	30" (75cm)	Violet	VETD-030R36	12 pcs	\$38.00
		0	30" (75cm)	Violet	VETD-030R30	12 pcs	\$38.00
X-LARGE FSL / C-16 30mm, 3/8 Circle		0	30" (75cm)	Violet	VETD-030R30	12 pcs	\$38.00

Needle Size	Needle Shape	Thread Size (USP)	Thread Length	Thread Color	Item #	Qty.	Price/box
<b>VET PDO</b> <span style="float: right;">Comparable to PDS II® * and Maxon® **</span>							
<b>LARGE</b> FS-1 / C-14 / C7 24mm, 3/8 Circle		2/0	30" (75cm)	Violet	VETD-230R24	12 pcs	\$38.00
		3/0	30" (75cm)	Violet	VETD-330R24	12 pcs	\$38.00
<b>MEDIUM</b> FS-2 / C-13 / C6 19mm, 3/8 Circle		4/0	30" (75cm)	Violet	VETD-430R19	12 pcs	\$38.00
<b>VET PGCL</b> <span style="float: right;">Comparable to Monocryl® * and Biosyn®**</span>							
<b>X-LARGE</b> CP-1 / GS-11 36mm, 1/2 Circle		0	30" (75cm)	Violet	VETQ-030R36-C	12 pcs	\$38.00
		0	30" (75cm)	Violet	VETQ-030R30-C	12 pcs	\$38.00
<b>LARGE</b> FS-1 / C-14 / C7 24mm, 3/8 Circle		2/0	30" (75cm)	Violet	VETQ-230R24-C	12 pcs	\$38.00
		3/0	30" (75cm)	Violet	VETQ-330R24-C	12 pcs	\$38.00
<b>MEDIUM</b> FS-2 / C-13 / C6 19mm, 3/8 Circle		4/0	30" (75cm)	Violet	VETQ-430R19-C	12 pcs	\$38.00
<b>VET Chromic Gut</b>							
<b>X-LARGE</b> CP-1 / GS-11 36mm, 1/2 Circle		0	30" (75cm)	Brown	VETC-030R36	12 pcs	\$32.00
		0	30" (75cm)	Brown	VETC-030R30	12 pcs	\$32.00
<b>LARGE</b> FS-1 / C-14 / C7 24mm, 3/8 Circle		2/0	30" (75cm)	Brown	VETC-230R24	12 pcs	\$32.00
		3/0	30" (75cm)	Brown	VETC-330R24	12 pcs	\$32.00
<b>MEDIUM</b> FS-2 / C-13 / C6 19mm, 3/8 Circle		4/0	30" (75cm)	Brown	VETC-430R19	12 pcs	\$32.00

**Topical Skin Adhesives**

GLUSEAL® Topical Skin Adhesives are liquid adhesive bandages intended to cover cuts, scrapes, burns and minor irritations of the skin and protect them from infection. Each 5ml multi-use kit contains enough tissue adhesive for up to 50 applications. Cost-effective for busy surgery practices. Available in multi-use and single-use versions. Made in Canada.



Item # **GLU-550**  
**Packing:** One Multi-Use Kit (1 x 5ml skin adhesive, 50 x micro-pipettes, 1 x autoclavable tray, 25 x foam tipped applicators)  
**Price:** \$130.00 /kit  
 \$510.00 /case of 4 kits



Item # **GLU-55012**  
**Packing:** Single-use kit  
 12 x 0.2ml /box  
**Price:** \$130 /box  
 \$510 /case of 4 boxes

**Disposable Skin Staplers**

EASYTAP® Disposable Skin Staplers are lever action staplers with an ergonomic design and strong construction to provide reliable performance during surgery. Each stapler is preloaded with 35 wide (6.5mm x 4.7mm) staples. Individually sealed and sterilized by EO Gas. Made in Korea.



Item # **SM-35W**  
**Packing:** 6 sealed staplers per box  
**Price:** \$79.50 /box  
 \$312.00 /case of 4 boxes

**Suturing Practice Skin Pads**

AD SURGICAL Suturing Practice Skin Pads are made from a durable high-grade silicone material that can be sutured repeatedly for practice. Each ultra-realistic pad mimics closely the 3 main layers of the skin (skin/ fat/ muscle) in texture and appearance. Available with or without pre-made wounds. Made in USA.

**5" x 7" pad with wounds (no case)**



Item # **T57-3SWP**  
**Packing:** 1 PC  
**Price:** \$49.00 each  
 \$475.00 /case of 10

**5" x 5" pad without wound (with case)**



Item # **T57-Skin**  
**Packing:** 1 PC  
**Price:** \$59.00 each  
 \$570.00 /case of 10

### Surgical Skin Markers (Violet Ink)


AD SURGICAL Surgical Skin Marker is designed to mark surgical incisions and anatomical sites on a patient's skin prior to surgical procedure to promote correct site surgery. Each sterile and ready-to-use surgical skin marker includes a 6" ruler and 6 blank labels. Two tip sizes available for different user and/or surgical procedure needs. Natural gentian violet plant extract ink is hypoallergenic and non-toxic. For single use only. Sterilized by gamma radiation.



Size	Item#	BOX OF 50
0.5mm Fine Tip	SKM-0005	\$49.50
1.0mm Regular Tip	SKM-1000	\$49.50
0.5mm + 1.0mm Tips	SKM-0510	\$49.50

### Aesthetic Skin Markers (Erasable White Ink)

AD SURGICAL Aesthetic Skin Marker has erasable white ink that can be easily wiped off without washing or scrubbing. Ideal for consultation appointments and/or non-surgical aesthetic procedures like Botox®, dermal fillers, Kybella® injection sites, and ultrasonic liposuction. Available only in 1.0mm (Regular) tip. For single use only. Sterilized by gamma radiation.



Size	Item#	BOX OF 50
1.0mm Tip	SKM-WH1000	\$49.50

### Cassette Sutures (for VET use only)

UNIFY® PDO Cassette Sutures are high-quality suture thread materials provided in a convenient and economical delivery system. All cassette sutures are intended for veterinary use only. Sterilized by EO Gas.



PDO				
0	15 meters	Violet	VCD-L0-15	\$54.00
2/0	15 meters	Violet	VCD-20-15	\$54.00
3/0	15 meters	Violet	VCD-30-15	\$54.00

### Suture Removal Scissors

AD SURGICAL Suture Removal Scissors are ideal for suture cutting and removal. The curved, blunt-tipped blades help locate and cut buried sutures without traumatizing adjacent sensitive tissue. Each scissors is 4 3/4" long and is fully autoclavable with a lifetime no-rust guarantee. Made in USA.



Item # SCR-B115  
 Packing: 1 pair  
 Price: \$45.00 each  
 \$252.00 /case of 6

### Winged IV Catheters (without safety)

SHARD Premium+ Winged IV Catheters (without safety) are designed to perform in delicate vasculature, delivering maximum efficiency and patient comfort. The smooth catheter-to-needle transition results in atraumatic insertion for the patient. Catheter features a flashback window for quick visualization of blood. Hub is color-coded for rapid identification of catheter size, and the winged design offers enhanced control and stability. Sterilized by EO Gas.



Size:	Item#	Box of 50:
18G x 1.25" (32mm)	#CAT-1832	\$59.50 /box
20G x 1" (25mm)	#CAT-2025	\$59.50 /box
20G x 1.25" (32mm)	#CAT-2032	\$59.50 /box
22G x 1" (25mm)	#CAT-2225	\$59.50 /box
24G x 0.75" (19mm)	#CAT-2419	\$59.50 /box

### Prefilled Saline IV Flush Syringes

SAFE-T-FLUSH Prefilled Normal Saline Flush Syringe (10ml) contains 0.9% Sodium Chloride USP. Compatible with all Luer-lock needle-free IV access systems. Single-use prefilled syringes help reduce medication errors and contamination risks associated with multi-dose vials. Sterile and for single use only.



Item # PFSFS10  
 Packing: 30 Prefilled Syringes / box  
 Price: \$25.00 /box  
 \$240.00 /case of 10 boxes

### Closed Safety IV Catheters (w/ Dual Port)


MERICATH Closed Safety IV Catheters (w/ Dual Port) is a complete system with tubing and extension set for administering IV therapy safely. It combines four individual devices in one: IV catheter, extension set, stabilization device, and two ports with a needle-free connector. This catheter system enhances patient safety and comfort. Sterilized by EO Gas.



Size: 22G X 1"  
 Item # MCC2250  
 Packing: 50 sets / box  
 Price: \$240.00 /box

### Blood Collection Sets (with Holder)

KBMED Safety Blood Collection Set with Holder features wide, flexible wings that provide greater control while positioning the needle. The safety design helps prevent needlestick injuries using a safety shield that pushes over the needle following venipuncture. Pre-attached holder meets OSHA single-use holder standard. Sterilized by EO gas.



Size:	Item#	Box of 25:
21G x 3/4"	#SBCN21-1H	\$35.00 /box
23G x 3/4"	#SBCN23-1H	\$35.00 /box
25G x 3/4"	#SBCN25-1H	\$35.00 /box

### Blood Collection Sets (w/o Holder)

KB MED Blood Collection Sets (without holder) allow users to collect blood draw samples with minimal risk due to their safety mechanism. They feature wide, flexible wings that provide greater control while positioning the needle, making one-handed application easier. The safety design helps prevent needlestick injuries using a safety shield that pushes over the needle following venipuncture. Sterilized by EO gas.



Size:	Item#	Box of 50:
21G x 3/4"	#SBCN21-1	\$45.00 /box
23G x 3/4"	#SBCN23-1	\$45.00 /box
25G x 3/4"	#SBCN25-1	\$45.00 /box

### Blood Collection Tube Holders

Blood Collection Tube Holders are made of a clear plastic for easy viewing of blood collection while providing a secure connection to blood collection sets, needles, and tubes when performing a direct draw into a blood collection tube. Designed to improve safety and compliance during routine blood collection procedures. Non-sterile and for single use only.



Item # **BCH-250**  
 Packing: 250 pcs /bag  
 Price: **\$25.00 /bag**  
**\$96.00 /case of 1000**

### JBP Ultra-thin Wall Nanoneedles

Ultra-thin Wall Nanoneedles by JBP are designed to reduce patient discomfort, minimize hematomas and limited swelling. Using a patented process called Micron Sharpness, the nanoneedles allow injection of higher viscosity solutions through a smaller entry to reduce pain and bleeding. Goal is to minimize the surface area so injectors have minimal impact on soft tissue and more satisfied patients. Sterilized by EO gas.

#### 30G



Size: 30G (4mm)  
 Item # **N-UTW3004**  
 Packing: 100 sterile capsules per box  
 Price: **\$75.00 /box**  
**\$730.00 /case of 10 boxes**

#### 33G



Size: 33G (13mm)  
 Item # **N-UTW3313**  
 Packing: 100 sterile capsules per box  
 Price: **\$75.00 /box**  
**\$730.00 /case of 10 boxes**

#### 34G



Size: 34G (8mm)  
 Item # **N-UTW3408**  
 Packing: 100 sterile capsules per box  
 Price: **\$85.00 /box**  
**\$830.00 /case of 10 boxes**

Size: 34G (4mm)  
 Item # **N-UTW3404**  
 Packing: 100 sterile capsules per box  
 Price: **\$85.00 /box**  
**\$830.00 /case of 10 boxes**

### JBP Ultra-thin Wall Nano Cannulas

Ultra-thin Wall Nano Cannulas by JBP are designed to reduce patient discomfort, minimize hematomas and limit swelling. Flexible cannula with blunt tip slides smoothly under the skin without causing trauma. Ideal for mesotherapy, placenta injection, Botox injection, collagen injection and local anesthetics. Each cannula comes packaged with a sharp punching needle. Sterilized by EO Gas.



Size:	Item#	Box of 24:
22G x 2" (50mm)	#CANN-2250	\$99.50 /box
22G x 2 3/4" (70mm)	#CANN-2270	\$99.50 /box
25G x 1 1/2" (40mm)	#CANN-2540	\$99.50 /box
25G x 2" (50mm)	#CANN-2550	\$99.50 /box
25G x 2 3/4" (70mm)	#CANN-2570	\$99.50 /box
27G x 1 1/2" (40mm)	#CANN-2740	\$99.50 /box
27G x 2" (50mm)	#CANN-2750	\$99.50 /box
30G x 1" (25mm)	#CANN-3025	\$99.50 /box

### Microsurgical Blades

SHARD® Microsurgical Blades are developed for minimally invasive procedures and microsurgery. These finer blades are suitable for use in many surgical disciplines. Made of high-grade carbon or stainless steel with an ultra-sharp cutting edge for clean incisions. High surface finishing ensures smooth texture. Each single-use blade is individually wrapped and sterilized by gamma radiation.



#61 (Chisel Tip)

#A600-M61  
12 blades / box  
**\$19.50 /box**



#62 (Chisel Tip)

#A600-M62  
12 blades / box  
**\$19.50 /box**



#64 (Round Tip)

#A600-M64  
12 blades / box  
**\$19.50 /box**



#65 (Lance Tip)

#A600-M65  
12 blades / box  
**\$19.50 /box**



#67 (Curved Tip)

#A600-M67  
12 blades / box  
**\$19.50 /box**



#69 (Sharp All Around)

#A600-M69  
12 blades / box  
**\$19.50 /box**



#90 (Spear Point Tip)  
*Ultra-Fine*

#A600-M90  
25 blades / box  
**\$29.50 /box**



#91 (Spear Point Tip)  
*Ultra-Fine*

#A600-M91  
25 blades / box  
**\$29.50 /box**

### Microsurgical Blade Handle

SHARD® Microsurgical Blade Handles feature a 4.5" (115mm) long round handle with serrated thumbs for better control. Made of surgical-grade stainless steel that is fully autoclavable. Compatible with all SHARD microsurgical blades.



#A600-MBH  
1 Piece  
**\$26.00 /piece**

### Premium Surgical Blades

SHARD® Premium+ Surgical Blades are available in 5 popular shapes and are used throughout all surgical disciplines including dermatology, oral surgery, orthopedics and reconstructive surgery. Computer-controlled grinding ensures a consistently sharp edge. Sterilized by gamma radiation.



#10 ITEM# **A600-P10-0**



#11 ITEM# **A600-P11-0**



#12 ITEM# **A600-P12-0**



#15 ITEM# **A600-P15-0**



#15C ITEM# **A600-P15-C**

Packing: 100 individually sealed blades per box  
Price: **\$34.00 /box**  
**\$320.00 /case of 10 boxes**

### Disposable Surgical Scalpels

Our economical SHARD® Disposable Scalpels incorporate a razor-sharp SS blade on an ergonomic plastic handle for excellent performance and control.



#10 ITEM# **A303-R10-0**



#11 ITEM# **A303-R11-0**



#12 ITEM# **A303-R12-0**



#15 ITEM# **A303-R15-0**



#15C ITEM# **A303-R15-C**



#20 ITEM# **A303-R20-0**

Packing: 100 individually sealed blades per box  
Price: **\$6.50 /box**  
**\$550.00 /case of 100 boxes**

### Disposable Dermal Curettes

SHARD® Premium+ Disposable Dermal Curettes are used for scraping skin lesions and growths such as warts or melanoma. It is a single-ended instrument available with a round loop tip. The stainless steel loop tip is sharp on the inside and blunt on the outside so the lesion can be removed, but surrounding skin is not damaged. Each curette is individually sealed and sterilized by gamma radiation.



Size 2mm	#A700-DC2	Size 5mm	#A700-DC5
Size 3mm	#A700-DC3	Size 7mm	#A700-DC7
Size 4mm	#A700-DC4		

Packing: 10 pcs/box  
Price: **\$29.50 /box**  
**\$280.00 /case of 10 boxes**

### Disposable Biopsy Punch

SHARD® Premium+ Disposable Biopsy Punch is ergonomically designed for routine dermal procedures. The razor-sharp seamless, stainless steel blade ensures safe, clean and easy sample extraction. Features a unique, easy-grip plastic handle with pressure assistance platform for superb control and comfort. Hollow handle allows easy specimen retrieval. Each punch is individually sealed and sterilized by gamma radiation.



Size 1.5mm	#A750-BP15	Size 4.0mm	#A750-BP40
Size 2.0mm	#A750-BP20	Size 5 mm	#A750-BP50
Size 2.5mm	#A750-BP25	Size 6 mm	#A750-BP60
Size 3.0mm	#A750-BP30	Size 8 mm	#A750-BP80
Size 3.5mm	#A750-BP35		

Packing: 10 pcs/box  
Price: **\$26.00 /box**  
**\$245.00 /case of 10 boxes**

### Disposable Dermal Biopsy Blade

SHARD® Premium+ Dermal Biopsy Blades are sterile, flexible scalpels designed for cutaneous surgery, including shave biopsies, saucerization of flat lesions, and leveling of pedunculated lesions. Our unique dermal blade design features a comfortable and protective finger guard along with an ultra-sharp, flexible blade. This flexibility allows the blade to be positioned at the correct angle for the intended procedure.



Item # **A750-DBB**  
Packing: 50 pcs/box  
Price: **\$105.00 /box**  
**\$400.00 /case of 4 boxes**

### Absorbable Hemostatic Dressings

BENACEL® Absorbable Hemostatic Dressing is a specially formulated gauze made of biocompatible oxidized cellulose and contains no harmful animal by-products or chemical additives. When applied to a wound, it produces hemostasis almost instantly by forming a gel-like substance upon contact with blood. BENACEL® is used to control bleeding during surgical procedures, or as a topical dressing on superficial open wounds. Made in USA.



# B01-0108	(0.5cm x 0.7cm plug)	8PCS/Box	<b>\$42.00 / Box</b>
# B01-0308	(0.6cm x 0.8cm plug)	8PCS/Box	<b>\$46.00 / Box</b>
# B01-0212	0.6"x0.6" (1.5cm x 1.5cm)	12PCS/Box	<b>\$55.00 / Box</b>
# B01-0510	2"x2" (5cm x 5cm)	10PCS/Box	<b>\$66.00 / Box</b>
# B02-1010	2"x4" (5cm x 10cm)	10PCS/Box	<b>\$148.00 / Box</b>

### Silicone Scar Removal Sheets

WOUND FREE® Silicone Scar Removal Sheets contain medical-grade silicone that helps hydrate scar tissue and mimics the natural barrier function of healthy skin to make scar softer, flatter and less visible. Each self-adherent sheet can last for up to 7-10 days of use. Color: Tan.



#SCS-T1R	1.5" x 2.4" (4cm x 6cm)	10 PCS / Box	<b>\$16.00 / box</b>
#SCS-T5R	1.2" x 6" (3cm x 15cm)	10 PCS / Box	<b>\$29.50 / box</b>

### Collagen Wound Dressing Pads

Made of 100% bovine collagen, WOUND FREE® Resorbable Collagen Wound Dressing Pads encourage the body to heal wounds quickly and effectively. Once applied to the wound, the collagen dressing forms a soft and comfortable gel on contact with the fluid that promotes healing throughout the day. The dressing can be precisely rolled or trimmed to fit the size of the wound. For medical and dental use. Sterile and for single use only.



# CWD-0505	2"x 2" (5cm x 5cm)	10 PCS / Box	<b>\$75.00</b>
		case of 5 Boxes	<b>\$360.00</b>

### Silicone Foam Dressings

WOUND FREE® Sterile Silicone Foam Dressings combine the superior absorption of foam and the gentle adhesion of silicone. The foam pad absorbs large amount of exudates. The silicone adhesive border provides gentle fixation, allowing the dressing to be removed or reapply when needed.



# FMD-S3	3"x3" (7.5cm x 7.5cm)	10PCS/Box	<b>\$18.00 / box</b>
# FMD-S4	4"x4" (10cm x 10cm)	10PCS/Box	<b>\$22.00 / box</b>
# FMD-S6	6"x6" (15cm x 15cm)	10PCS/Box	<b>\$36.00 / box</b>

### Adhesive Island Dressings

WOUND FREE® Sterile Adhesive Island Dressing is a gentle soft cloth bordered gauze dressing indicated for wounds with light to moderate exudation. One-step application is an excellent alternative to tape and gauze to help save user time. The non-adherent central gauze pad provides absorption of light to moderate drainage.



# ISD-P0708	2.75" x 3.15" (7cm x 8cm)	55PCS/Box	<b>\$12.00 / box</b>
# ISD-P1010	4" x 4" (10cm x 10cm)	40PCS/Box	<b>\$14.00 / box</b>

### Super Absorbent Dressings

WOUND FREE® Sterile Super Absorbent Dressing is a new generation of advanced wound dressing. It can be used as either a primary or secondary dressing with the purpose of using its extremely absorbent core for exudate management. The self-adherent super absorbent polymer dressing core effectively soaks up high amounts of wound exudate from pressure and venous ulcers, surgical incisions, burns, donor grafts and chronic wounds.



# SAD-AS4	4"x4" (10cm x 10cm)	10PCS/Box	<b>\$12.00 / box</b>
-----------	---------------------	-----------	----------------------

### Post-op Therapy Gel Packs

AD SURGICAL Hot/Cold Therapy Gel Packs are frequently used by cosmetic and plastic surgery patients to reduce swelling and bruising. They are also used by patients after oral surgeries including wisdom teeth extractions, reconstructive jaw surgery, and implant surgery. Our reusable Hot/Cold Gel Packs stay pliable when frozen and remain cold for over an hour.



<b>Cold Gel Pack</b>	#G400-CLR	10 PCS / Pack	<b>\$23.00</b>
<b>3" x 5" Clear</b>		50 PCS / Case	<b>\$110.00</b>
<b>Hot + Cold Gel Pack</b>	#GB425-BL	12 PCS / Pack	<b>\$32.00</b>
<b>3" x 6" Galaxy Blue</b>		48 PCS / Case	<b>\$125.00</b>

### Reusable Cold Gel Packs

AD SURGICAL Reusable Cold Gel Pack helps relieve pain and swelling after a surgical procedure or during the initial stage of an acute injury. 4" diameter with soft cloth backing provides comfortable cold therapy for up to 20 minutes. Cold gel pack is commonly used for dermal fillers, chemical peels, facial injections, oral surgery and other surgical procedures.



#T430-BU	4" Round	10 PCS / Bag	<b>\$25.00 / Bag</b>
#T430-CS-BU	4" Round	50 PCS / Bag	<b>\$120.00 / Bag</b>

### Soft Stretch Facial Wraps

AD SURGICAL Post-Op Soft Stretch Facial Wraps by COOL JAW® feature an exclusive, hand-free patented design. Wraps are made with a spandex fabric that is extremely soft on the patient's skin. They feature a continuous outer sleeve that extends cheekbone to cheekbone to ensure full coverage and compression over the surgical area. This sleeve is specifically designed to fit our gel packs and allows patients to easily insert (and replace) gel packs without having to remove the wrap. The wrap features two rows of Velcro closure at the top of the head to provide a universal, adjustable fit.



Face Wrap With Chin Cup, Black  
Item# W800C-BLK  
Price: \$16.00 /piece  
\$186.00 /12 pieces



Face Wrap Without Chin Cup, White  
Item# W800-WHT  
Price: \$16.00 /piece  
\$186.00 /12 pieces

### Universal Surgical Drapes

AD SURGICAL Sterile Universal Surgical Drapes are made of an absorbent and completely impervious PE/ PP material to prevent strike-through and subsequent contamination of the surgical site. Each drape has a peel-able adhesive tape and can be used alone or in combination, giving flexibility to individual procedural needs. Individually wrapped and sterilized by EO Gas. Latex free.



#### HALF DRAPE WITH TAPE

Sterile Half Drape with Adhesive Tape

- 60in x 44in ( 152cm x 112cm )
- Absorbent, Impervious PE/PP Material

Item # A400-HFD  
Packing: 10 individually sealed packs per box  
Price: \$34.00 /box  
\$190.00 /case of 6 boxes

#### MINOR PROCEDURE DRAPE with adhesive aperture

Sterile Drape with Adhesive Tape

- 22in x 25in ( 56cm x 64cm )
- Aperture Size: 2.5in ( 6cm ) Diameter
- Absorbent, Impervious PE/PP Material

Item # A400-APD  
Packing: 20 individually sealed packs per box  
Price: \$30.00/box  
\$345.00 /case of 12 boxes

#### UTILITY DRAPE WITH TAPE

Sterile Drape with Adhesive Tape

- 22in x 25in ( 56cm x 64cm )
- Absorbent, Impervious PE/PP Material

Item # A400-UTD  
Packing: 20 individually sealed packs per box  
Price: \$30.00/box  
\$345.00 /case of 12 boxes

### Nitrile Surgical Gloves (Non-Latex, Powder-Free)

SUPERFIT® Sterile Nitrile Surgical Gloves are specially designed for people with sensitive skin. These hypoallergenic gloves are made from synthetic nitrile rubber that are non-latex and powder-free. The well-designed curved finger and anatomic shape help reduce hand fatigue, and the textured surface in the grip area provides for maximum gripping.




Available Sizes	Item # for Box
Size 6.0	A560-N60
Size 6.5	A560-N65
Size 7.0	A560-N70
Size 7.5	A560-N75
Size 8.0	A560-N80
Size 8.5	A560-N85

Packing: 50 pairs per box  
Price: \$65.00 /box  
\$378.00 /case of 6 boxes

### Sterile Hose & Cable Covers

EZ DRAPE® Sterile Hose and Cable Covers are user friendly and effective. Simply insert the handpiece or other surgical instrument through the ring and pull the plastic sleeve down. Sterile hose cover helps maintain a clean and safe surgical environment. Individually wrapped and sterilized by EO Gas. Latex free.

Ring diameter: 2.2" Length: 60" long



Item # A400-1000  
Packing: 24 individually sealed packs per box  
Price: \$64.00/box  
\$364.00 /case of 6 boxes

### Latex Surgical Gloves (Powder-Free)

SUPERFIT® Sterile Latex Surgical Gloves are made from ergonomic molds that help reduce hand fatigue during operation. Powder-free reduces the risk of latex allergy. Straight finger design and wet donning ability make it very easy to don either with dry or damp hands. Each pair of Left/ Right fitted gloves are sealed in a sterile pouch and sterilized by gamma / ebeam radiation.



Available Sizes	Item # for Box
Size 6.0	A520-LTX60
Size 6.5	A520-LTX65
Size 7.0	A520-LTX70
Size 7.5	A520-LTX75
Size 8.0	A520-LTX80

Packing: 50 pairs per box  
Price: \$49.50 /box  
\$288.00 /case of 6 boxes

## Sterile Oral Surgery Drape Packs

O-PAK® Sterile Oral Surgery Drape Pack helps maintain a clean and safe surgical environment to prevent postoperative infections and cross contamination. Ideal for dental implant surgery or other aseptic procedures. Simplify procedural setup with 14 essential disposable products included in one convenient package. Sterilized by EO Gas. Latex free.

Each pack contains all the items listed below :  
( Items are not sold separately )

- 2 Procedure gowns, universal size
- 2 Facemasks with ties, 2 Bouffant caps
- Absorbent cotton scrub towels, 20 x 40 cm,  
1 Surgical instrument tray
- 1 Impervious patient drape with tape,  
112 x 152 cm
- 2 Table or mayo stand covers with adhesive tape,  
64 x 43 cm
- 1 Absorbent O/R towel
- 1 Saliva ejector, 1 Surgical aspirator tip, 2 Cotton tip  
applicators, 25 Gauze swabs
- 4 Light handle/ chair button barrier films
- 1 Medical wrapping sheet, 90 x 90 cm



Item # **A200-3000**

Packing: 6 individually sealed packs per box

Price: **\$120.00 /box**

**\$472.00 /case of 4 boxes**

## Oral Tissue Adhesives

PERIACRYL® Oral Tissue Adhesive is a biocompatible and non-toxic cyanoacrylate tissue adhesive that is useful in securing sutures, stabilizing membranes, covering bone grafting material post-extraction, or dressing donor and recipient sites. Violet coloring allows for better visibility. Available in multi-use and single-use versions. Made in Canada.




**Multi-Use Kit  
(High Viscosity)**

Item# **GLU-PERI50**

Packing: 5 ml bottle  
50 pipettes + 1 Tray

Price: **\$130.00 /box**

**\$512.00 /case of 4**



**Single-Use 12 Pack  
(High Viscosity)**

Item# **GLU-PERI12**

Packing: 12 x (0.2 ml) / box

Price: **\$130.00 /box**

**\$512.00 /case of 4**

## Thermoplastic Surgical Guides for Pilot-Drilling

EZ STENT® Thermoplastic Surgical Guide for Pilot-Drilling is our original patented design that allows fabrication of an accurate surgical guide in minutes using hot water alone. The template has a titanium drill sleeve with inner diameter 2.3mm that is compatible with the pilot drill of most dental implant systems. EZ STENT® should be cold sterilized prior to usage and is for single use only.

### SINGLE TOOTH SET

EZ STENT® Single Tooth Set allows accurate preparation of an osseous site in a tooth-bound region. Titanium drill sleeve with inner Ø2.3mm is compatible with the pilot drill of most dental implant systems.



Item # **D015-003-P**

Contents:

- 5 Thermoplastic templates for Single-Tooth sites  
- Sleeve height 10mm, sleeve inner Ø2.3mm
- 5 Guide pins ( stainless steel ), Length 20mm, Ø2.3mm
- 1 Drill for dental laboratory ( steel ), Ø2.3mm

Price: **\$160.00 /set**

**\$560.00 /case of 4 sets**

### FREE END SET

EZ STENT® Free End Set allows accurate preparation of an osseous site in the most distal region. Titanium drill sleeve with inner Ø2.3mm is compatible with the pilot drill of most dental implant systems.



Item # **D015-004-P**

Contents:

- 5 Thermoplastic templates for Free-End Situations  
- Sleeve height 10mm, sleeve inner Ø2.3mm
- 5 Guide pins ( stainless steel ), Length 20mm, Ø2.3mm
- 1 Drill for dental laboratory ( steel ), Ø2.3mm

Price: **\$160.00 /set**

**\$560.00 /case of 4 sets**

### 2-UNIT BRIDGE SET

EZ STENT® 2-Unit Bridge Set allows the surgeon to accurately locate and drill the holes of two consecutive implants with optimal spacing, angulation and parallelism. Titanium drill sleeves with inner Ø2.3mm are compatible with the pilot drill of most dental implant systems.



Item # **D015-008-2U**

Contents:

- 4 Thermoplastic templates ( titanium/ polymer )  
- 2 Tooth-Bound and 2 Free End Templates  
- Sleeve height 6mm, sleeve inner Ø2.3mm
- 8 Guide pins ( stainless steel ), Length 20mm, Ø2.3mm
- 1 Drill for dental laboratory ( steel ), Ø2.3mm

Price: **\$160.00 /set**

**\$560.00 /case of 4 sets**

## Pilot Surgical Drills (for EZ Stent®)

Our Ø2.3mm Pilot Surgical Drills are perfect fit for the EZ STENT® Pilot-Drill Sets. Feature wide-band, easy to see, laser-etched depth markings at 8mm, 10mm, 11.5mm, 13mm and 16mm. Made in USA.



Item # **D015-023L**

Content:

- 1 Long pilot surgical drill  
( stainless steel )  
Length: 38mm

Price: **\$65.00 /drill**



Item # **D015-023S**

Content:

- 1 Short pilot surgical drill  
( stainless steel )  
Length: 30mm

Price: **\$65.00 /drill**



### Universal Implant Drill Stopper Set

EZ STOPPER® Universal Implant Drill Stopper Set is a set of medical grade, autoclavable silicone tubings. Each tubing can be accurately measured and cut to the desired length in seconds using the built-in lockable metal stop in the autoclavable storage tray. The trimmed silicone stoppers will stretch-fit over most surgical drills to provide precise drill depth control during dental implant placement. Each silicone stopper is color-coded to fit the corresponding surgical drill sizes as follows:

- RED Ø 1.8 – 2.4mm
- GREEN Ø 2.5 – 3.4mm
- YELLOW Ø 3.5 – 4.4mm
- BLUE Ø 4.5 – 5.6mm



**Autoclavable**

EZ Stopper Starter Kit:

Item #D017-STPS

Packing: 1 set (includes 4 color coded tubing)

Price: **\$160.00 /set**

Silicone Tubing Refills:

Red ( Item #D017-STP18 )

Green ( Item #D017-STP25 )

Yellow ( Item #D017-STP35 )

Blue ( Item #D017-STP45 )

Packing: 4 PCS per refill pack

Price: **\$28.00 /pack**



### Implant Abutment PTFE Fillers

i-PLUG® implant Abutment PTFE Filler is a semi-soft and packable filler material that can be used to seal an abutment hole. It can also be easily retrieved by our Filling & Removal (FR) Instrument when needed. The one-size PTFE rod (Ø1.8mm x 50mm) can be cut to custom fit the abutment hole of any screw-retained or cement-retained crown. i-PLUG® is fully autoclavable to prevent cross-contamination and to minimize wastage of any unused material.



**Autoclavable**

Item # D030-050

Packing: 1 box (of 15)

Price: **\$49.50 /box**

**\$192.00 /case of 4 boxes**

### Filling & Removal (FR) Instrument

Double-sided instrument with both a plugger end and a barbed retrieval end designed for handling abudment fillers. Made of a surgical grade stainless steel and is fully autoclavable and rust resistance.



Item # D030-017

Packing: 1 double-sided instrument

Price: **\$39.50**

### AUGMA Degranulation Kit

AUGMA Degranulation Kit is intended for removing granulation tissue before augmentation. Ideal for removing all the remnants of soft tissue which are well attached to the bone and grab them out without harming the bone. The debridement process is done quickly and effectively. Use the bur without irrigation while removing soft tissue and use it with irrigation when grinding bones. Multi-task drills can be used in implantology, periodontal surgery, endodontic surgery, and root removal. The speed can be adjusted according to the kind of the procedure performed hence preventing emphysema. Each kit includes one bur each of the following sizes: 1mm, 2.5mm, 3mm, 3.5mm.



**Autoclavable**



Item # D035-D-KIT

Packing: 1 Kit

Price: **\$230.00 /Kit**

### Implant Impression Trays

IMPLANT TRAY® Implant Impression Tray has an innovative patented design that features a wide array of clear and easily removable panels to accommodate open-tray impression techniques. This eliminates guesswork and minimizes time needed to drill out holes in the impression tray. Made in USA.



**Large Upper**

Item# D020-U3L

**Medium Upper**

Item# D020-U2M

**Small Upper**

Item# D020-U1S

**Large Lower**

Item# D020-L3L

**Medium Lower**

Item# D020-L2M

**Small Lower**

Item# D020-L1S

Packing: 12 trays per box

Price: **\$30.00 /box**

#### Assortment Starter Package

Item # D020-COMBO

Contents: 6 of each sizes ( Large Upper, Medium Upper, Small Upper, Large Lower, Medium Lower, Small Lower )

Packing: 36 trays per box

Price: **\$84.00 /box**

### Synthetic Bioactive Bone Grafts

UNIGRAFT® Bioactive Bone Graft is made of a synthetic bioactive glass material. It offers an osteogenic framework for new bone growth and bolsters a higher level of osteoblastic activity to enhance bone regeneration. UNIGRAFT® is completely resorbable and transforms into natural bone with time. Made in USA.



- Bone Graft (200-600um)** Packing: 5 vials per box
- 0.5 gram Item# U03-0105 Price: **\$182.00 /Box**
  - 1.0 gram Item# U03-0205 Price: **\$292.00 /Box**
  - 2.0 gram Item# U03-0305 Price: **\$396.00 /Box**

### Mineralized Allograft In Syringes

OSSEOSEAL™ bone in syringe is a mixture of allograft mineralized cortical and cancellous powder in a convenient syringe with a filter cap! Curved syringe for better access to site. The cancellous and cortical blend forms an osteoconductive scaffold providing volume enhancement and effective site development for successful ridge preservation and dental implant placement. Sterile and for single use only.



- Bone Graft (250-800um)**
- 0.3cc Item# OS-03cc Price: **\$44.00 /Pack**
  - 0.5cc Item# OS-05cc Price: **\$66.00 /Pack**
  - 1.0cc ( 2 of 0.5cc syringes in 1 Pack)  
Item# OS-10cc Price: **\$100.00 /Pack**

### Bovine Bone Grafts

OSSEOSEAL™ Bovine Bone Graft is a xenograft produced by removal of organic components from bovine bone. It provides a supportive structure for osteoconduction. The porous bone mineral matrix has macro and microscopic structures that mimic human bone. The formation and ingrowth of new bone at implantation site is favored, due to its trabecular architecture, interconnecting macro and micro pores.



- Powder (200-1000um)** Packing: 1 vial per box
- 0.5 gram Item# OSX-05CC Price: **\$38.50 /Box**
  - 1.0 gram Item# OSX-10CC Price: **\$55.00 /Box**
  - 4.0 gram Item# OSX-40CC Price: **\$132.00 /Box**

### BOND APATITE Bone Graft Cements

BOND APATITE by Augma Materials is a composite graft, made of Biphasic Calcium Sulfate mixed with HA granules in a controlled particle size distribution, intended to fill or augment a broad range of bony defects. It fully resorbs and is replaced by bone during the healing process. It allows minimal invasive surgical procedure, reduced treatment time and easy manipulation of the graft material for the clinician.



- Item# HA-BB-PACK  
Packing: 1cc Syringe / Box  
Price: **\$135.00 /Box**  
**\$1080.00 /Case of 8**

### Mineralized Allograft in Vials

REVOSS™ Mineralized Human Allograft in Vials is a sterile processed human allograft tissue ideal for a variety of surgical applications. It comprises 50/50 mineralized cortical / cancellous bone powder. All tissue was recovered, processed, stored and distributed for use in accordance with the standards of AATB and FDA requirements.



- 0.25 CC Item# RB5-1025 Price: **\$44.00 /Box**
- 0.5 CC Item# RB5-1050 Price: **\$75.00 /Box**
- 1.0 CC Item# RB5-1100 Price: **\$90.00 /Box**
- 2.0 CC Item# RB5-1200 Price: **\$150.00 /Box**
- 2.5 CC Item# RB5-1250 Price: **\$166.00 /Box**

### Demineralized Allograft Putty

REVOSS™ Demineralized Allograft Putty is a demineralized bone matrix (DBM) allograft made of donated human tissue. It is packaged in an easy-to-use prefilled syringe for direct application with no need for pre-treatment in a variety of surgical applications. All tissue was recovered, processed, stored and distributed for use in accordance with the standards of AATB and FDA requirements.



#### Demineralized Allograft Putty in Syringes

- 0.7 cc Item# RP2-0070 Price: **\$136.00 /Piece**
- 1.0 cc Item# RP2-0100 Price: **\$158.00 /Piece**
- 2.5 cc Item# RP2-0250 Price: **\$370.00 /Piece**

### Resorbable Bone Substitutes

BICERA® Resorbable Bone Substitute is a bioceramic material that is composed of 60% hydroxyapatite (HAP) and 40% beta-tricalcium (β-TCP). It can be used in common dental procedures such as sinus lifting, periodontal bone fillings, and filling of peri-implant defects. BICERA provides a bone void filler that resorbs and is replaced by bone during the healing process.



#### Bicera HA/TCP Granule 1-2mm

- 2cc Item# BC-G22 Price: **\$160.00 /Box**
- 5cc Item# BC-G02 Price: **\$240.00 /box**
- 10cc Item# BC-G07 Price: **\$400.00 /Box**
- Block Item# BC-B01 Price: **\$210.00 /Box**

### Allograft Bone Fiber Plugs

REVOSS® Matrix Bone Fiber Plugs are composed of a human allograft mixture of demineralized cortical bone fiber and mineralized cancellous bone particulates. The unique, compressed cylinder shape of REVOSS Matrix allows for easy insertion and expansion within the bone void, eliminating the need for particulate mixing or membranes to facilitate bone regeneration.



#### 11mm x 7.5mm Bone Insert

- Item# RM1-0101  
Packing: 1 Piece / Box  
Price: **\$96.00 /Box**  
**\$470.00 /Case of 5**

## Attachable Intraoral Wound Dressings

ORA-AID® is an adhesive intraoral patch used to protect a wound such as an extraction socket, tissue grafting & implant site, and ulcers. It is easy to use, just peel off and stick on the wound site. The non-resorbable polymer helps cover the affected area to protect it until soft tissue closes and healing begins.



### 25mm x 15mm Intraoral Dressing

Item# AG-202A

Packing: 20 Pieces / Box

Price: **\$99.00 /Box**

**\$485.00 /Case of 5 boxes**

## Resorbable Collagen Plugs

HEALIAID® Resorbable Collagen Plugs are designed to deliver dependable wound dressing at extraction and biopsy sites. Shaped for ease of use and compressible to suit open wounds, these resorbable dressings are inexpensive and easy to work with. When used to dress fresh extraction sockets, the collagen plugs protect exposed surfaces from external irritants and further injury. The Type-1 purified bovine collagen is non-pyrogenic, and typically reabsorbed by the body within 2 to 4 weeks after placement. Each collagen plug comes sterile and is ready for use without preparation.



### 10mm x 20mm Collagen Plug

Item# SA001020

Packing: 10 Plugs / Box

Price: **\$105.00 /Box**

**\$515.00 /Case of 5 boxes**

## Bilayer Collagen Matrix

HEALIAID® Bilayer Collagen Matrix is designed for soft tissue regeneration and sealing extraction sockets in the oral cavity. It is composed of two layers: the dense outer layer provides stability while allowing open healing; the spongy inner layer supports blood clot stabilization and promotes early vascularization and good soft-tissue ingrowth.

The resorbable dressing is indicated for gaining keratinized tissue and for recession coverage. It is also used for sealing extraction sockets with preserved buccal walls, which helps retain bone grafting materials. The collagen matrix provides an alternative to autologous soft-tissue grafts. Painful harvesting of tissue is avoided, benefiting both patients and clinicians. The matrix is 3mm thickness.



### 10mm x 10mm (for sealing extraction sockets)

Item# HADB101030

Packing: 5 Pieces / Box

Price: **\$220.00 /Box**

### 15mm x 20mm (for sealing molar extraction sockets)

Item# HADB152030

Packing: 1 Piece / Box

Price: **\$65.00 /Box**

**\$300.00 /Case of 5**

### 25mm x 30mm (for soft tissue regeneration)

Item# HADB253030

Packing: 1 Piece / Box

Price: **\$140.00 /Box**

**\$680.00 /Case of 5**

### 30mm x 40mm (for soft tissue regeneration)

Item# HADB304030

Packing: 1 Piece / Box

Price: **\$200.00 /Box**

**\$980.00 /Case of 5**

## Resorbable PLA / PGA Barrier Membranes

CYTOFLEX RESORB® Resorbable PLA / PGA Barrier Membrane consists of three integral layers designed to promote gingival tissue attachment and prevent fibroblast down growth. Made of biocompatible poly (lactide-co-glycolide) copolymers, CYTOFLEX RESORB® is non-pyrogenic, non-immunogenic and maintains a barrier framework for over two months after implantation. The resorbable membrane does not require pre-soaking before administration, and completely dissolves within six months. Made in USA.

### 12mm x 24mm

#### Resorbable Membrane

Item# C03-0301

Packing: 1 per pack

Price: **\$60.00 /pack**

**\$294.00 /case of 5**



### 20mm x 25mm

#### Resorbable Membrane

Item# C03-0201

Packing: 1 per pack

Price: **\$88.00 /pack**

**\$432.00 /case of 5**



### 30mm x 40mm

#### Resorbable Membrane

Item# C03-0101

Packing: 1 per pack

Price: **\$124.00 /pack**

**\$608.00 /case of 5**



## Resorbable Collagen Membranes

COLLAFLEX™ Resorbable Collagen Membrane is derived from purified porcine pericardium. The non-crosslinked membrane retains the natural microfibril structure of the pericardium, providing excellent biocompatibility and extended barrier function for guided bone regeneration applications. Highly drapable and can be easily fixated with sutures and bone tacks.

### 15mm x 20mm

#### Collagen Membrane

Item# C07-0101

Packing: 1 per pack

Price: **\$75.00 /pack**

**\$368.00 /case of 5**



### 20mm x 30mm

#### Collagen Membrane

Item# C07-0201

Packing: 1 per pack

Price: **\$100.00 /pack**

**\$490.00 /case of 5**



### 30mm x 40mm

#### Collagen Membrane

Item# C07-0301

Packing: 1 per pack

Price: **\$166.00 /pack**

**\$815.00 /case of 5**



### ePTFE Barrier Membranes

CYTOFLEX TEF-GUARD® ePTFE Barrier Membrane is a non-resorbable membrane made of microporous ePTFE material. It is designed to resist the penetration of fibroblasts and bacteria, while simultaneously allowing the exchange of interstitial fluids through the membrane. This enhances gingival tissue attachment and provides a favorable protected environment for neo-vascularization and re-population of osteoblasts in the bony defects. Made in USA.



**12mm x 24mm  
ePTFE Membrane**

Item# C04-0101  
Packing: 1 per pack  
Price: **\$54.00 /pack**  
**\$240.00 /case of 5**



**25mm x 30mm  
ePTFE Membrane**


Item# C04-0201  
Packing: 1 per pack  
Price: **\$72.00 /pack**  
**\$306.00 /case of 5**

### Titanium-Reinforced ePTFE Barrier Membranes

CYTOFLEX Ti-ENFORCED® Titanium-reinforced ePTFE Barrier Membrane is a multi-layer, non-resorbable membrane. It is a barrier material with a titanium frame enclosed within two layers of ePTFE material. The titanium reinforcement is intended for space and shape maintenance to contain bone grafting material and to minimize graft movement during wound healing. The Ti-Enforced™ membranes are textured on both surfaces to improve handling and adhesion between tissue and membrane. Made in USA.



**12mm x 24mm**  
Item# C05-0301  
Packing: 1 per pack  
Price: **\$109.00 /pack**  
**\$535.00 /case of 5**



**19mm x 26mm**  
Item# C05-0901  
Packing: 1 per pack  
Price: **\$148.00 /pack**  
**\$726.00 /case of 5**




**14mm x 24mm**  
Item# C05-0501  
Packing: 1 per pack  
Price: **\$136.00 /pack**  
**\$668.00 /case of 5**



**20mm x 25mm**  
Item# C05-1301  
Packing: 1 per pack  
Price: **\$180.00 /pack**  
**\$882.00 /case of 5**




**25mm x 30mm**  
Item# C05-1701  
Packing: 1 per pack  
Price: **\$200.00 /pack**  
**\$980.00 /case of 5**



**32mm x 40mm**  
Item# C05-4101  
Packing: 1 per pack  
Price: **\$245.00 /pack**  
**\$1200.00 /case of 5**

### Titanium Mesh Membranes

CYTOFLEX® Titanium Mesh Membrane is made of a pure, implantable titanium mesh with precision profile. It is designed for use to ensure three-dimensional reconstruction of alveolar bone defects and to facilitate bone augmentation through adequate fixation of the augmentation material. Made in USA.



**12mm x 25mm  
Mesh Non-Resorbable  
Membrane**

Item# C02-4101  
Packing: 1 per pack  
Price: **\$48.00 /pack**  
**\$236.00 /case of 5**



**20mm x 25mm  
Mesh Non-Resorbable  
Membrane**

Item# C02-4301  
Packing: 1 per pack  
Price: **\$60.00 /pack**  
**\$294.00 /case of 5**



**25mm x 30mm  
Mesh Non-Resorbable  
Membrane**

Item# C02-4401  
Packing: 1 per pack  
Price: **\$60.00 /pack**  
**\$294.00 /case of 5**

### Sterile Oral Exam Packs



STERILE i-PAK® is an economical 3-in-1 intraoral exam pack that eliminates cross-contamination, improves efficiency and saves on sterilization costs. Each i-pak® includes 1 disposable dental mirror with anti-fog coating and 1 patented disposable double-ended Explorer & Probe. Probe is marked 3-6-9 mm universal style, tapered with round tip for excellent readability and patient comfort. Safe for use on titanium implants. The explorer is made with surgical grade stainless steel with an ultra-sharp tip. Good for routine exam, emergency treatment, off-site visit and high volume oral exam performed by schools, government and other agencies. Sterilized by EO Gas. Latex free.

Item # D002-009-P  
Packing: 50 individually sterilized packs per box  
Price: **\$86.00 /box**  
**\$1320.00 /case of 16 boxes**

**AD Surgical**  
1296 Kifer Road, Suite 608,  
Sunnyvale CA 94086 USA

**FREE SAMPLES**  
for most products!



**AD Surgical**

Shop online: [www.AD-SURGICAL.com](http://www.AD-SURGICAL.com)

Order by phone: 888.841.8481

## A Simple Promise

Our promise is to deliver the best customer experience by offering:

- Same day shipping on most orders
- A 90-day risk free return policy on all products
- A team of friendly and knowledgeable customer support specialists to ensure your complete satisfaction

**We look forward to serving you soon.**



Scan to learn more about us

*Please keep this catalog for reference.*

All prices are in U.S. dollars and are subject to applicable taxes according to state and local laws. Availability, prices and delivery rates are subject to change without notice. Please visit our website for the most updated pricing & product details.