

UNIFY® Suture Brand Name Cross Reference Chart

	UNIFY®	ETHICON®	MEDTRONIC®
ABSORBABLE	PGA	VICRYL®	POLYSORB®
	PGCL	MONOCRYL®	BIOSYN®
	RAPID PGA	VICRYL® RAPIDE	VELOSORB® FAST
	PDO	PDS® / PDS® II	MAXON®
	CHROMIC GUT	CHROMIC GUT	CHROMIC GUT
	PLAIN GUT	PLAIN GUT	PLAIN GUT
	BARBED PDO	STRATAFIX® PDS	V-LOC® 180
NON-ABSORBABLE	NYLON	ETHILON®	MONOSOF® / DERMALON®
	POLYPROPYLENE	PROLENE®	SURGIPRO®
	PTFE	N/A	N/A
	SILK	PERMA-HAND® SILK	SOFSILK®
	POLYESTER	ETHIBOND® EXCEL	TI-CRON®
	STAINLESS STEEL	STAINLESS STEEL	STEEL

@2023 AD Surgical. All rights reserved. UNIFY and UNIFY logo are trademark of AD Surgical. Other third party brands are trademarks of their respective owners.