

















UNIFY® Suture Needle Comparison Chart


▼ Reverse Cutting Needles













ETHICON/ UNIFY	MEDTRONIC	SHARPOINT / LOOK	
CP	GS-12	C14	 40mm, 1/2 Circle
CP-1	GS-11	-	 36mm, 1/2 Circle
PSL	P-17	-	 30mm, 3/8 Circle
FSL	C-16	FSL	 30mm, 3/8 Circle
PS-1	P-14	PC-33	 24mm, 3/8 Circle
FS-1	C-14	C7	 24mm, 3/8 Circle
PS-2	P-12	PC-31	 19mm, 3/8 Circle
FS-2	C-13	C6	 19mm, 3/8 Circle
PS-3	P-11	PC34	 16mm, 3/8 Circle
FS-3	C-12	C22	 16mm, 3/8 Circle
P-3	P-13	C3	 13mm, 3/8 Circle
P-1	P-10	C1	 11mm, 3/8 Circle




Taper Point Needles






ETHICON/ UNIFY	MEDTRONIC	SHARPOINT / LOOK	
CTX	GS-25	CTX	 48mm, 1/2 Circle
CT	GS-24	T10	 40mm, 1/2 Circle
CT-1	GS-21	T8	 36mm, 1/2 Circle
SH	V-20	C21	 26mm, 1/2 Circle
RB-1	CV-23	T28	 17mm, 1/2 Circle
TF	CV-22	T22	 13mm, 1/2 Circle
MICRO BV-130-5	-	-	 6.4mm, 1/2 Circle
MICRO BV-130-4	MICRO MV-100-4	-	 5mm, 1/2 Circle
MICRO BV-100-3	MICRO MV-70-3	-	 4mm, 1/2 Circle




Conventional Cutting Needles




ETHICON/ UNIFY	MEDTRONIC	SHARPOINT / LOOK	
PC-3	PC-11	PC-3	 16mm, 3/8 Circle
PC-1	PC-13	-	 13mm, 3/8 Circle




Straight Cutting Needles






ETHICON/ UNIFY	MEDTRONIC	SHARPOINT / LOOK	
KS	SC-2	KS	 60mm



Micropoint Spatula Needles



ETHICON/ UNIFY	MEDTRONIC	SHARPOINT / LOOK	
SP-13	-	-	 13mm, 3/8 Circle
TG140-8 (Single)	SE-140-8	-	 6.4mm, 3/8 Circle
TG140-8 (Double)	SE-140-6	-	 6.4mm, 3/8 Circle